

Urgent Field Safety Notice

ACHC-20-02.A.OUS.DM

February 2020

Dimension® clinical chemistry system

Eltrombopag Interference with Dimension® Total Bilirubin (TBI) Flex® reagent cartridge

Our records indicate that your facility may have received the following product:

Table 1: Dimension® Chemistry Products affected:

Assay	Test Code	Catalog Number	Siemens Material Number (SMN)	Lot Number
Total Bilirubin	TBI	DF167	10444957	ALL

Reason for Correction

The purpose of this communication is to inform you of an interference identified with the products indicated in Table 1 above and provide instructions on actions that your laboratory must take.

Siemens Healthcare Diagnostics has become aware that the United Kingdom Medicines and Healthcare products Regulatory Agency published an alert to healthcare professionals informing them that laboratory tests for bilirubin should be monitored for patients who take the drug eltrombopag due to the potential for discordant results. Eltrombopag may be used in the treatment of thrombocytopenia and/or aplastic anemia. Siemens spiking studies have shown a positive bias in Total Bilirubin (TBI) results at a therapeutic concentration of eltrombopag. Interference has not been observed with the Direct Bilirubin (DBI) assay.

Table 2 below reflects eltrombopag Interference with Dimension® Total Bilirubin (TBI) assay based on Siemens internal testing.

Table 2. Eltrombopag Interference with Dimension Total Bilirubin Assay

Analyte	Analyte Concentration mg/dL [µmol/L]	Eltrombopag Concentration µg/mL [µmol/L]	Bias (%)
TBI	0.8 [13.7]	25 [56.5]	89.3
TBI	22 [376]	25 [56.5]	3.5

The "Limitations of the Procedure" section in the Instructions For Use (IFU) for the Dimension TBI assay will be updated to indicate: *Use of this assay is not recommended for patients undergoing treatment with eltrombopag due to the potential for falsely elevated results.*

Information related to eltrombopag provided in this letter supersedes the information in the current Dimension TBI IFU until the IFU is updated. Siemens will communicate once the IFU has been updated.

Risk to Health

For patients taking eltrombopag, the potential exists for the misinterpretation of total bilirubin levels, which may confound investigations for the etiology of hyperbilirubinemia. Potential clinical impact would be mitigated by correlation to clinical symptomology and additional laboratory testing including other markers of liver function (e.g. alanine aminotransferase, aspartate aminotransferase, alkaline phosphatase, and/or lactate dehydrogenase). Siemens is not recommending a review of previously generated results.

Actions to be Taken by the Customer:

- For patients on eltrombopag therapy, use of Dimension TBI is not recommended.
- Please review this letter with your Medical Director.
- Complete and return the Field Correction Effectiveness Check Form attached to this letter within 30 days.
- If you have received any complaints of illness or adverse events associated with the products listed in Table 1, immediately contact your local Siemens Customer Care Center or your local Siemens technical support representative.

Please retain this letter with your laboratory records and forward this letter to those who may have received this product.

We apologize for the inconvenience this situation may cause. If you have any questions, please contact your Siemens Customer Care Center or your local Siemens Technical Support representative.

Dimension is a trademark of Siemens Healthcare Diagnostics.

FIELD CORRECTION EFFECTIVENESS CHECK

Eltrombopag Interference with the Dimension® Total Bilirubin (TBI) Assay.

This response form is to confirm receipt of the enclosed Siemens Healthcare Diagnostics Urgent Field Safety Notice ACHC-20-02.A.OUS.DM dated February 2020 regarding eltrombopag Interference with the Dimension® Total Bilirubin (TBI) Flex® reagent cartridge. Please read the question and indicate the appropriate answer.

Return this completed form to Siemens Healthcare Diagnostics as per the instructions provided at the bottom of this page.

1. I have read and understood the Urgent Field Safety Notice instructions provided in this letter. Yes No

Name of person completing questionnaire: _____

Date: _____

Title: _____

Institution: _____ Instrument Serial Number: _____

Street: _____

City: _____ State: _____

Phone: _____ Country: _____

Please send a scanned copy of the completed form via email to your local Siemens Healthineers technical support representative by email at CruinnFSNGroup@cruinn.ie or by fax to 01-6297401

If you have any questions, contact your local Siemens Healthineers technical support representative.