



---

## **Urgent Field Safety Notice**

*Recall*

**concerning**

Anti-IA2 ELISA (IgG), order no. EA 1023-9601 G, lot no. E220215BA

---

3 June 2022

**From:**

EUROIMMUN Medizinische Labordiagnostika AG  
Seekamp 31  
23560 Lübeck  
Germany  
www.euroimmun.com

**To:**

Users and distributors

**Identification of the medical/IVD products concerned:**

Anti-IA2 ELISA (IgG), order no. EA 1023-9601 G, lot no. E220215BA

Dear Customer,

EUROIMMUN has initiated a field corrective action for the product Anti-IA2 ELISA (IgG), order no. EA 1023-9601 G, lot no. E220215BA. This notification contains important information for your immediate attention.

**Description of the problem and determined cause:**

Due to customer complaints in relation to invalid positive controls due to low optical density (OD) values for calibrator 3 during processing of the Anti-IA2 ELISA (IgG) (order no. EA 1023-9601 G) of lot E220215BA, we have initiated an investigation to confirm the existence of the problem, identify its cause and determine an appropriate corrective measure.

So far, the investigation has revealed that the calibrators 3 included in the aforementioned product lot might contain a concentration of IA2 antibodies that is too low, leading to a significantly reduced measurement signal (OD value) for calibrator 3 compared to the target value indicated on the quality control certificate valid for this lot.

The investigation also showed that this problem does not affect the other calibrators (1, 2, 4, 5, 6) and the controls included in the relevant test kit lot.

As a consequence of this significant measurement signal deviation for calibrator 3, the positive control contained in the test kit gives an invalid result. An invalid positive control indicates that the respective test processing and the results obtained for the analysed patient samples are invalid. It is possible that diagnoses are delayed due to these invalid test runs, but there is no risk of incorrect diagnosis. If the calibration is valid, the positive control shows as valid and the test results can be used. Previous analyses of patient samples in which the affected lot was used and in which that was the case must therefore not be reanalysed.



**Measures to be taken:**

The fault encountered for calibrator 3 has no influence on the qualitative results of the test. In the quantitative evaluation using the 450 nm calibration curve in accordance with the instructions for use, the faulty calibrator 3 has no negative influence on the quantitative results obtained for samples in the negative or weak positive concentration range. However, it will cause invalid results for the positive control included in the test kit and incorrect quantitative results for patient samples with a higher concentration if no evaluation method with automatic curve correction is used for the 405 nm calibration curve.

Please make sure that all remaining and potentially faulty calibrator 3 of the indicated lot are no longer used in your laboratory. We kindly ask that you fill in the template on page 3 to confirm that you have received this safety notice and that you send it by fax to the following number: +49 (0) 451 2032 7065 immediately, **latest until the 13th of June 2022.**

EUROIMMUN will replace the calibrator 3 free of charge. Please contact our colleague Anke Fletemeyer in the complaints department by email to [a.fletemeyer@euroimmun.de](mailto:a.fletemeyer@euroimmun.de).

**Information to be passed on:**

This notice must be forwarded to all users and distributors of the above-mentioned product.

Thank you for your cooperation! We apologise for any inconvenience this may cause.

For further information, please do not hesitate to contact any of the following contact persons at EUROIMMUN.

**Contact persons:**

Product Management Endocrinology  
Fax: +49 (0) 451 2032 7065  
E-Mail: [endocrinology-pm@euroimmun.de](mailto:endocrinology-pm@euroimmun.de)

PRRC-V Immunobiochemical Tests  
Dr. Christian Krüger  
Tel.: +49 (0) 151 22617145  
Fax: +49 (0) 451 2032 100  
E-Mail: [c.krueger@euroimmun.de](mailto:c.krueger@euroimmun.de)

---

*Signature / Person Responsible for Regulatory Compliance- Vigilance (PRRC-V)*

EUROIMMUN Medizinische Labordiagnostika AG  
Seekamp 31  
23560 Lübeck  
Germany



---

**Please send back the distributor reply form as specified on the document!**



## Field Safety Notice

### Distributor Reply Form

Mandatory fields are marked with \*

1. Field Safety Notice (FSN) information	
FSN Reference number*	2022-4
FSN Date*	03.06.2022
Product/ Device name*	Anti-IA2 ELISA (IgG)
Product Code(s)	EA 1023-9601 G
Batch/Serial Number (s)	E220215BA

2. Distributor Details	
Company Name*	
Account Number	
Address*	
Shipping address if different to above	
Contact Name*	
Title or Function	
Telephone number*	
Email*	

3. Return acknowledgement to Sender	
Email	<a href="mailto:endocrinology-pm@euroimmun.de">endocrinology-pm@euroimmun.de</a>
Distributor Helpline	Product recall 2022-4
Postal Address	Seekamp 31, 23560 Lübeck
Web Portal	<a href="http://www.euroimmun.de">www.euroimmun.de</a>
Deadline for returning the Distributor/Importer reply form*	13.06.2022

4. Distributors (Tick all that apply)		
<input type="checkbox"/>	*I confirm the receipt, the reading and understanding of the Field Safety Notice.	Distributor to complete or enter N/A
<input type="checkbox"/>	I have checked my stock and quarantined inventory	Distributor to enter quantity and date
<input type="checkbox"/>	I have identified customers that received or may have received this device	
<input type="checkbox"/>	I have attached customer list	
<input type="checkbox"/>	I have informed the identified customers of this FSN	Date of communication:
<input type="checkbox"/>	I have received confirmation of reply from all identified customers	
<input type="checkbox"/>	I have destroyed affected reagents (calibrator 3) - enter number destroyed and date complete.	Add quantity, Lot/Serial Number/Date of Destruction
<input type="checkbox"/>	Neither I nor any of my customers has any affected devices in inventory	
Print Name*		
Signature*		
Date *		



---

Please send the completed Field Safety Notice Form immediately to:

EUROIMMUN Medizinische Labordiagnostika AG  
Seekamp 31  
23560 Lübeck  
Germany  
Fax: +49 451 2032 7065  
Email: [endocrinology-pm@euroimmun.de](mailto:endocrinology-pm@euroimmun.de)