03 juuli 2023

**KIIRELOOMULINE OHUTUSTEATIS – MMS-23-4760**

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Tootja seerianumber (SRN): CH-MF-000026539

Kataloogi nr (VIITENR):	Tootekirjeldus	UDI-DI
80203UNS01-035	Alaris™ Gateway Workstation 0x3x5 Configuration without intelligent base unit	N/A – Released prior to UDI being assigned.
80203UNS01-235	Alaris™ Gateway Workstation 2x3x5 Configuration without intelligent base unit	N/A – Released prior to UDI being assigned.
80203UNS01-30	Alaris™ Gateway Workstation 3x0 Configuration without intelligent base unit	N/A – Released prior to UDI being assigned.
80203UNS01-32	Alaris™ Gateway Workstation 3x2 Configuration without intelligent base unit	N/A – Released prior to UDI being assigned.
80203UNS01-33	Alaris™ Gateway Workstation 3x3 Configuration without intelligent base unit	N/A – Released prior to UDI being assigned.
80203UNS01-34	Alaris™ Gateway Workstation 3x4 Configuration without intelligent base unit	N/A – Released prior to UDI being assigned.
80203UNS01-50	Alaris™ Gateway Workstation 5x0 Configuration without intelligent base unit	N/A – Released prior to UDI being assigned.
80203UNS01-52	Alaris™ Gateway Workstation 5x2 Configuration without intelligent base unit	N/A – Released prior to UDI being assigned.
80203UNS01-53	Alaris™ Gateway Workstation 5x3 Configuration without intelligent base unit	N/A – Released prior to UDI being assigned.
80203UNS01-54	Alaris™ Gateway Workstation 5x4 Configuration without intelligent base unit	N/A – Released prior to UDI being assigned.
80203UNS01-70	Alaris™ Gateway Workstation 7x0 Configuration without intelligent base unit	N/A – Released prior to UDI being assigned.
80203UNS01-72	Alaris™ Gateway Workstation 7x2 Configuration without intelligent base unit	N/A – Released prior to UDI being assigned.
80203UNS01-73	Alaris™ Gateway Workstation 7x3 Configuration without intelligent base unit	N/A – Released prior to UDI being assigned.
80203UNS01-74	Alaris™ Gateway Workstation 7x4 Configuration without intelligent base unit	N/A – Released prior to UDI being assigned.
80203UNS01-92	Alaris™ Gateway Workstation 9x2 Configuration	N/A – Released prior to UDI being assigned.
80203UNS02-035	Alaris™ Gateway Workstation 0x3x5 Configuration with the wired Ethernet communication capability and bar code reader interface	N/A – Released prior to UDI being assigned.
80203UNS02-235	Alaris™ Gateway Workstation 2x3x5 Configuration with the wired Ethernet communication capability and bar code reader interface	N/A – Released prior to UDI being assigned.
80203UNS02-30	Alaris™ Gateway Workstation 3x0 Configuration with the wired Ethernet communication capability and bar code reader interface	N/A – Released prior to UDI being assigned.
80203UNS02-32	Alaris™ Gateway Workstation 3x2 Configuration with the wired Ethernet communication capability and bar code reader interface	N/A – Released prior to UDI being assigned.
80203UNS02-33	Alaris™ Gateway Workstation 3x3 Configuration with the wired Ethernet communication capability and bar code reader interface	N/A – Released prior to UDI being assigned.
80203UNS02-34	Alaris™ Gateway Workstation 3x4 Configuration with the wired Ethernet communication capability and bar code reader interface	N/A – Released prior to UDI being assigned.
80203UNS02-50	Alaris™ Gateway Workstation 5x0 Configuration with the wired Ethernet communication capability and bar code reader interface	N/A – Released prior to UDI being assigned.
80203UNS02-52	Alaris™ Gateway Workstation 5x2 Configuration with the wired Ethernet communication capability and bar code reader interface	N/A – Released prior to UDI being assigned.
80203UNS02-53	Alaris™ Gateway Workstation 5x3 Configuration with the wired Ethernet communication capability and bar code reader interface	N/A – Released prior to UDI being assigned.
80203UNS02-54	Alaris™ Gateway Workstation 5x4 Configuration with the wired Ethernet communication capability and bar code reader interface	N/A – Released prior to UDI being assigned.
80203UNS02-70	Alaris™ Gateway Workstation 7x0 Configuration with the wired Ethernet communication capability and bar code reader interface	N/A – Released prior to UDI being assigned.
80203UNS02-72	Alaris™ Gateway Workstation 7x2 Configuration with the wired Ethernet communication capability and bar code reader interface	N/A – Released prior to UDI being assigned.

80203UNS02-73	Alaris™ Gateway Workstation 7x3 Configuration with the wired Ethernet communication capability and bar code reader interface	N/A – Released prior to UDI being assigned.
80203UNS02-74	Alaris™ Gateway Workstation 7x4 Configuration with the wired Ethernet communication capability and bar code reader interface	N/A – Released prior to UDI being assigned.
80203UNS02-92	Alaris™ Gateway Workstation 9x2 Configuration with the wired Ethernet communication capability and bar code reader interface	N/A – Released prior to UDI being assigned.
80223UNS02-30	Alaris™ Gateway Workstation 3x0 Configuration with wired Ethernet connectivity to ACE	10885403478659
80223UNS02-035	Alaris™ Gateway Workstation 0x3x5 Configuration with wired Ethernet connectivity to ACE	10885403478789
80223UNS02-32	Alaris™ Gateway Workstation 3x2 Configuration with wired Ethernet connectivity to ACE	10885403478666
80223UNS02-33	Alaris™ Gateway Workstation 3x3 Configuration with wired Ethernet connectivity to ACE	10885403478673
80223UNS02-34	Alaris™ Gateway Workstation 3x4 Configuration with wired Ethernet connectivity to ACE	10885403478680
80223UNS02-50	Alaris™ Gateway Workstation 5x0 Configuration with wired Ethernet connectivity to ACE	10885403478697
80223UNS02-52	Alaris™ Gateway Workstation 5x2 Configuration with wired Ethernet connectivity to ACE	10885403478703
80223UNS02-53	Alaris™ Gateway Workstation 5x3 Configuration with wired Ethernet connectivity to ACE	10885403478710
80223UNS02-54	Alaris™ Gateway Workstation 5x4 Configuration with wired Ethernet connectivity to ACE	10885403478727
80223UNS02-70	Alaris™ Gateway Workstation 7x0 Configuration with wired Ethernet connectivity to ACE	10885403478734
80223UNS02-72	Alaris™ Gateway Workstation 7x2 Configuration with wired Ethernet connectivity to ACE	10885403478741
80223UNS02-73	Alaris™ Gateway Workstation 7x3 Configuration with wired Ethernet connectivity to ACE	10885403478758
80223UNS02-74	Alaris™ Gateway Workstation 7x4 Configuration with wired Ethernet connectivity to ACE	10885403478765
80223UNS02-92	Alaris™ Gateway Workstation 9x2 Configuration with wired Ethernet connectivity to ACE	10885403478772
80223UNS02-235	Alaris™ Gateway Workstation 2x3x5 Configuration with wired Ethernet connectivity to ACE	10885403478796
80203UNS03-30	Alaris™ Gateway Workstation 3x0 Configuration with a wireless communications interface	N/A – Released prior to UDI being assigned.
80203UNS03-32	Alaris™ Gateway Workstation 3x2 Configuration with a wireless communications interface	N/A – Released prior to UDI being assigned.
80203UNS03-33	Alaris™ Gateway Workstation 3x3 Configuration with a wireless communications interface	N/A – Released prior to UDI being assigned.
80203UNS03-34	Alaris™ Gateway Workstation 3x4 Configuration with a wireless communications interface	N/A – Released prior to UDI being assigned.
80203UNS03-50	Alaris™ Gateway Workstation 5x0 Configuration with a wireless communications interface	N/A – Released prior to UDI being assigned.
80203UNS03-52	Alaris™ Gateway Workstation 5x2 Configuration with a wireless communications interface	N/A – Released prior to UDI being assigned.
80203UNS03-53	Alaris™ Gateway Workstation 5x3 Configuration with a wireless communications interface	N/A – Released prior to UDI being assigned.
80203UNS03-54	Alaris™ Gateway Workstation 5x4 Configuration with a wireless communications interface	N/A – Released prior to UDI being assigned.
80203UNS03-70	Alaris™ Gateway Workstation 7x0 Configuration with a wireless communications interface	N/A – Released prior to UDI being assigned
80203UNS03-72	Alaris™ Gateway Workstation 7x2 Configuration with a wireless communications interface	N/A – Released prior to UDI being assigned
80203UNS03-73	Alaris™ Gateway Workstation 7x3 Configuration with a wireless communications interface	N/A – Released prior to UDI being assigned
80203UNS03-74	Alaris™ Gateway Workstation 7x4 Configuration with a wireless communications interface	N/A – Released prior to UDI being assigned
80300UNS01-035	Alaris™ Gateway Workstation 0x3x5 Configuration without intelligent base unit	10885403478611
80300UNS01-235	Alaris™ Gateway Workstation 2x3x5 Configuration without intelligent base unit	10885403478635
80300UNS01-30	Alaris™ Gateway Workstation 3x0 Configuration without intelligent base unit	10885403478239

80300UNS01-32	Alaris™ Gateway Workstation 3x2 Configuration without intelligent base unit	10885403478246
80300UNS01-33	Alaris™ Gateway Workstation 3x3 Configuration without intelligent base unit	10885403478253
80300UNS01-34	Alaris™ Gateway Workstation 3x4 Configuration without intelligent base unit	10885403478260
80300UNS01-50	Alaris™ Gateway Workstation 5x0 Configuration without intelligent base unit	10885403478277
80300UNS01-52	Alaris™ Gateway Workstation 5x2 Configuration without intelligent base unit	10885403478284
80300UNS01-53	Alaris™ Gateway Workstation 5x3 Configuration without intelligent base unit	10885403478291
80300UNS01-54	Alaris™ Gateway Workstation 5x4 Configuration without intelligent base unit	10885403478307
80300UNS01-70	Alaris™ Gateway Workstation 7x0 Configuration without intelligent base unit	10885403478314
80300UNS01-72	Alaris™ Gateway Workstation 7x2 Configuration without intelligent base unit	10885403478321
80300UNS01-73	Alaris™ Gateway Workstation 7x3 Configuration without intelligent base unit	10885403478338
80300UNS01-74	Alaris™ Gateway Workstation 7x4 Configuration without intelligent base unit	10885403478345
80300UNS01-92	Alaris™ Gateway Workstation 9x2 Configuration without intelligent base unit	10885403478352
80300UNS02-035	Alaris™ Gateway Workstation 0x3x5 Configuration with the wired Ethernet communication capability and bar code reader interface	10885403478628
80300UNS02-235	Alaris™ Gateway Workstation 2x3x5 Configuration with the wired Ethernet communication capability and bar code reader interface	10885403478642
80300UNS02-30	Alaris™ Gateway Workstation 3x0 Configuration with the wired Ethernet communication capability and bar code reader interface	10885403478369
80300UNS02-32	Alaris™ Gateway Workstation 3x2 Configuration with the wired Ethernet communication capability and bar code reader interface	10885403478376
80300UNS02-33	Alaris™ Gateway Workstation 3x3 Configuration with the wired Ethernet communication capability and bar code reader interface	10885403478383
80300UNS02-34	Alaris™ Gateway Workstation 3x4 Configuration with the wired Ethernet communication capability and bar code reader interface	10885403478390
80300UNS02-50	Alaris™ Gateway Workstation 5x0 Configuration with the wired Ethernet communication capability and bar code reader interface	10885403478406
80300UNS02-52	Alaris™ Gateway Workstation 5x2 Configuration with the wired Ethernet communication capability and bar code reader interface	10885403478413
80300UNS02-53	Alaris™ Gateway Workstation 5x3 Configuration with the wired Ethernet communication capability and bar code reader interface	10885403478420
80300UNS02-54	Alaris™ Gateway Workstation 5x4 Configuration with the wired Ethernet communication capability and bar code reader interface	10885403478437
80300UNS02-70	Alaris™ Gateway Workstation 7x0 Configuration with the wired Ethernet communication capability and bar code reader interface	10885403478444
80300UNS02-72	Alaris™ Gateway Workstation 7x2 Configuration with the wired Ethernet communication capability and bar code reader interface	10885403478451
80300UNS02-73	Alaris™ Gateway Workstation 7x3 Configuration with the wired Ethernet communication capability and bar code reader interface	10885403478468
80300UNS02-74	Alaris™ Gateway Workstation 7x4 Configuration with the wired Ethernet communication capability and bar code reader interface	10885403478475
80300UNS02-92	Alaris™ Gateway Workstation 9x2 Configuration with the wired Ethernet communication capability and bar code reader interface	10885403478482

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