

URGENT FIELD SAFETY NOTICE

Commercial Name: HandyVac™, Guedel Airways, Universal Connector, Exchange bags for UnoMeter™ 500 (see Attachment 2 for commercial name and full details)

Issue Date: 02nd October 2019
 Original Notice Revised Notice Revision No.: Rev. 1

Reference/REF No: See Attachment 2: Product List

FSCA Ref.: 2019-008-011MA & 013MA

Type of Action: Recall/Product Disposal

Please note that this action only applies to specific product codes of the above-mentioned products.

DESCRIPTION OF THE PROBLEM	<p>ConvaTec has voluntarily initiated a recall of specific product codes of the above-mentioned products.</p> <p>Internal assessment of product packaging integrity has confirmed that these devices are not meeting our expectations or those of our customers. Transportation testing conducted on the product in its current market packaging failed, confirming the potential for a breach in the sterile barrier.</p> <p>ConvaTec has not received any reports of incidents related to the packaging seal issue.</p> <p>Using a non-sterile device on a patient may expose the patient to infectious agents increasing the patient risk of developing infection.</p>
-----------------------------------	---

INFORMATION ON AFFECTED DEVICES	<p>HandyVac™ HandyVac™ is a single-use sterile low vacuum wound drainage system for post-operative surgical wound drainage, to drain wound fluids and prevent formation of hematomas, which can delay healing, cause pain or support infection. HandyVac™ consists of a drainage tube connected to a low-vacuum bellows and a drainage bag to collect the drained fluids. It is supplied as a standalone system or together with wound drains or wound drains that are pre-connected to a trocar.</p> <p>Guedel Airways Single-use sterile oropharyngeal airways to maintain a patent airway in unconscious or anaesthetized patients by keeping the airway passage free.</p> <p>Universal Connector Sterile connector with CH 04-18 Christmas Tree connector on one side and Luer slip port on the other side.</p> <p>Exchange Bags for UnoMeter™ 500 Sterile exchange bags for use with the UnoMeter™ 500 hourly diuresis monitoring system. The system is used to monitor urine output in critically ill patients. The affected exchange bags are used to collect larger amounts of urine over a prolonged period. Once the bag is full, it is usually replaced with a new bag if the system is not replaced as a whole.</p>
--	---

PRODUCT IDENTIFICATION PROCEDURE	<p><u>Confirmation of Specific Product Codes and Lot Numbers:</u></p> <p>This issue is limited to specific product codes as identified in Attachment 2.</p> <p>Only the identified product part codes within this notice may have a potential breach in the sterile barrier packaging.</p> <p>For this reason and to address any potential risk of harm, all the affected products should not be used.</p> <p>The only way to identify affected product is by comparing Product code/REF and LOT/Batch number to the recalled product list (see Attachment 1). There is no other discernable difference between affected and unaffected product.</p> <p>See Attachment 1 for examples of package labelling that highlights the location of the product code and LOT/Batch number on the device labels. These are located on the primary packaging and the shipping carton. The product code (reference number) is preceded by the word 'REF' and the LOT/Batch number is preceded by the symbol, LOT.</p>
---	--

DISTRIBUTOR ACTIONS

1	Verify if you have affected product in your inventory and identify which of your customers have received potentially affected product using Attachment 2: Product List .
2	<p>If you have distributed this product:</p> <ul style="list-style-type: none"> • To other Wholesalers: <ul style="list-style-type: none"> ○ Please forward this Field Safety Notice, including attachments, to them and ask that they follow these Distributor Actions. • To End Users <ul style="list-style-type: none"> ○ Please forward this Field Safety Notice, including attachments, to your end-users and ask that they follow the End User Actions Section (below).
3	<p>If you have received returned product from your customers that they have identified or if you have any stock on hand which you have verified as included in the scope of this Field Safety Notice (as described in Attachment 2: Product List).</p> <ul style="list-style-type: none"> • Immediately stop distributing and quarantine all affected products. • Contact ConvaTec Customer Service for further instructions on how to return affected product and to arrange credit. <p>Products will not be automatically replaced. For replacement product(s), please place new orders with ConvaTec using the usual process. ConvaTec will supply new products as soon as possible.</p>
4	<p>Please complete the Field Safety Notice <u>Distributor Response Form</u> and return to your Customer Service as soon as possible. Include details of product to be returned, if applicable.</p> <p>Please send a complete list of all consignees to your Customer Service as soon as possible. Please indicate if you have sent this Field Safety Notice, including attachments, to your consignees. This information is not required but important to allow ConvaTec to perform effectiveness checks.</p> <p>Please contact your Customer Service if you have any additional questions.</p>

PLEASE PROVIDE A COMPLETED RESPONSE AS SOON AS POSSIBLE

RETAILER ACTIONS

1	<p>For Retailers with a physical location (e.g. a “brick and mortar” store)</p> <ul style="list-style-type: none"> • Immediately stop distributing and quarantine all affected products • Please post the Field Safety Notice including Attachment 2 in a conspicuous location in your store. <p>For Retailers, which have a digital location (e.g. an eRetailer)</p> <ul style="list-style-type: none"> • Please forward Attachment 2: Product List showing affected Product Codes/REF No. and LOT Numbers via physical mail or e-mail to End-users.
2	<p>If you have received returned product from your customers that they have identified or if you have any stock on hand which you have verified as included in the scope of this Field Safety Notice (as described in Attachment 2: Product List, contact Customer Services Centre for further instructions on how to return affected product and to arrange credit.</p> <p>Products will not be automatically replaced. For replacement product (s), please place new orders with ConvaTec using the usual process. ConvaTec will supply new products as soon as possible.</p>
3	<p>Please complete the Field Safety Notice Distributor Response Form and return to your Customer Services as soon as possible. Include details of product to be returned, if applicable.</p> <p>Please send a complete list of all consignees to your Customer Service as soon as possible. Please indicate if you have sent this Field Safety Notice, including attachments, to your consignees.</p> <p>This information is not required but important to allow ConvaTec to perform effectiveness checks. Please contact your Customer Service if you have any additional questions.</p>

PLEASE PROVIDE A COMPLETED RESPONSE AS SOON AS POSSIBLE

END USER ACTIONS (CLINICIANS, HOSPITALS, PATIENTS, OTHERS)

1	<p>Our records show that you have taken delivery of affected product. Please stop the use of all affected devices as defined in this document.</p>
2	<p>This issue is limited to specific product codes and lots. Check stock and ensure that all affected devices that you have in stock are quarantined.</p> <p>If you do not know where the product code and/or lot number are located, please refer to Attachment 1: Examples of Labelling.</p> <p>To verify whether your lot is associated with the product included in the scope of this Field Safety Notice use Attachment 2: Product List showing affected Product Codes/REF No. and LOT Numbers.</p>
3	<p>Please complete the enclosed, Field Safety Notice ‘End User Response Form’ and forward it to your distributor as soon as possible.</p> <p>Include details of product to be returned, if applicable.</p> <p>This information is required to allow ConvaTec to perform effectiveness checks.</p>
4	<p>Contact your distributor (or ConvaTec Customer Service if you have purchased the product directly from ConvaTec) to arrange return of affected products, if applicable, and to arrange credit. For replacement products, please place new orders with your distributor or ConvaTec using the usual route and procedure.</p> <p>Please contact your Customer Service if you have any additional questions.</p>


PLEASE PROVIDE A COMPLETED RESPONSE AS SOON AS POSSIBLE

TRANSMISSION OF THIS FIELD SAFETY NOTICE	<p>This notice needs to be passed on to all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred. (As appropriate)</p> <p>Please transfer this notice to other organisations on which this action has an impact. (As appropriate)</p> <p>Please maintain awareness on this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action.</p> <p>Please report all device-related incidents to the manufacturer, distributor or local representative, and the national Competent Authority if appropriate, as this provides important feedback.</p>
---	---

ConvaTec is committed to providing quality products and services to our customers and we sincerely apologise for any inconvenience this notice may cause.

The relevant National Authorities have been advised about this Field Safety Corrective Action.

AUTHORISATION

<u>Name</u>	<u>Title</u>	<u>Address</u>	
Diane Baker	Director, CCC Franchise Regulatory Affairs and Quality Assurance	ConvaTec Limited, First Avenue, Deeside Industrial Park, Deeside, CH5 2NU, U.K.	
<u>Date</u>	02 October 2019	<u>Signature</u>	



Regional ConvaTec Customer Service Contact:

Austria, Belgium, Czech Republic, Estonia, France, Germany, Ireland, Israel, Netherlands, Poland, Spain, Switzerland & Turkey

Tel: +41 (0) 52 630 54 01

Email: ccc.customerservice@convatec.com

Denmark

Tel: +45 48 167 030

Email: customerservicenordic@convatec.com

Finland

Tel: +358 20 7659 630

Email: mail.fi@convatec.com

Italy

Tel: +39 06 89789928

Email: CustomerServiceItaly@Convatec.com

Norway

Tel: +47 22 686 095

Email: customerservicenordic@convatec.com

Sweden

Tel: +46 (0)42 332010

Email: customerservicenordic@convatec.com

United Kingdom

Tel: +44 (0) 1244 832206

Email: Unomedical-uk.customerservice@convatec.com



FIELD SAFETY NOTICE DISTRIBUTOR RESPONSE FORM
PLEASE COMPLETE AND RETURN by Fax/Email

Consignee of the device:

Consignee Account No:	
Consignee Name:	
Consignee Address:	

The following products, HandyVac™, Guedel Airways, Universal Connector and Exchange bags For UnoMeter™ 500 have been distributed to your facility:

Invoice #	Sales Order #	Product Code / REF No.	SAP Code	LOT No.	Quantity Delivered (pieces)

Distributors (Tick all that apply and give details, where applicable)		
<input type="checkbox"/>	I confirm the receipt, the reading and understanding of the Field Safety Notice.	
<input type="checkbox"/>	I have checked my stock and quarantined inventory	Add details to Table 1
<input type="checkbox"/>	I have identified customers that received or may have received this device	
<input type="checkbox"/>	I have attached customer list	Add details to Table 2
<input type="checkbox"/>	I have informed the identified customers of this FSN	Date sent:
<input type="checkbox"/>	I have received confirmation of reply from all identified customers	Attach responses
<input type="checkbox"/>	Neither I nor any of my customers has any affected devices in inventory	

It is important that your organisation takes the actions detailed in the FSN and confirms that you have received the FSN. Your organisation's reply is the evidence we need to monitor the progress of the corrective actions.



Table 1. Quarantined Inventory: Record quantity (pieces) for each LOT to be returned.

LOT No.	Units on Hand	LOT No.	Units on Hand	LOT No.	Units on Hand	LOT No.	Units on Hand

Table 2. Customer List: Please provide details of affected HandyVac™, Guedel Airways, Universal Connector and Exchange bags For UnoMeter™ 500 products that were distributed to your customers.

Customer Name	Product Code / REF No.	SAP Code	LOT No.	Quantity (pieces)

FORM Completed and Returned From:	
Name (CAPITAL LETTERS):	
Position:	
Company Name:	
Address:	
Phone No:	
Signature:	
Date (dd/mmm/yyyy):	



FIELD SAFETY NOTICE CUSTOMER RESPONSE FORM
PLEASE COMPLETE AND RETURN by Fax/Email

Consignee of the device:

Consignee Account No:	
Consignee Name:	
Consignee Address:	

The following products, HandyVac™, Guedel Airways, Universal Connector and Exchange bags For UnoMeter™ 500 have been distributed to your facility:

Invoice #	Sales Order #	Product Code / REF No.	SAP Code	LOT No.	Quantity Delivered (pieces)

Customer action undertaken on behalf of Healthcare Organisation (Tick all that apply)		
<input type="checkbox"/>	I confirm receipt of the Field Safety Notice and that I read and understand its content.	
<input type="checkbox"/>	I performed all actions requested by the FSN.	
<input type="checkbox"/>	The information and required actions have been brought to the attention of all relevant users and executed.	
<input type="checkbox"/>	I have checked my stock and quarantined inventory	Add details to Table 1
<input type="checkbox"/>	No affected devices are available for return	

It is important that your organisation takes the actions detailed in the FSN and confirms that you have received the FSN. Your organisation's reply is the evidence we need to monitor the progress of the corrective actions.

Table 1. Quarantined Inventory: Record quantity (pieces) for each LOT to be returned.

LOT No.	Units on Hand	LOT No.	Units on Hand	LOT No.	Units on Hand	LOT No.	Units on Hand











FORM Completed and Returned From:	
Name (CAPITAL LETTERS):	
Position:	
Company Name:	
Address:	
Phone No:	
Signature:	
Date (dd/mmm/yyyy):	

List of Attachments

Attachment 1: Examples of Labelling and Affected Product

Attachment 2: EU Product List showing affected Product Codes/REF No. and LOT Numbers

Attachment 1: Examples of Labelling and Affected Product

For the AFFECTED product	
<p>...the shipper labelling will include the product code and lot number.</p> <p>See example labels to the right.</p> <p>Product Code (REF/SKU)</p> <p>Lot Number</p>	<div style="border: 2px solid yellow; padding: 5px; display: inline-block; margin-bottom: 10px;">REF 84061181</div> <p> GB Unometer bag 2000 S Unometerpåse DK Unometerpose 2000 FIN Unometerpussi NO Unometerpose 2000 Made in Belarus </p> <p style="text-align: right; font-size: 24px;">2 X 100</p> <div style="display: flex; justify-content: space-between; align-items: center;"> <div style="text-align: center;">  </div> <div style="text-align: center;">  </div> <div style="font-size: 8px;"> ConvaTec Limited First Avenue Deeside Industrial Park Deeside, Flintshire, CH5 2NU, UK </div> </div> <div style="display: flex; justify-content: space-between; align-items: center; margin-top: 10px;"> <div style="text-align: center;">  2012-04  2017-03 </div> <div style="border: 2px solid blue; padding: 5px; display: inline-block;"> LOT 632574 </div> </div> <p style="text-align: center; font-size: 8px;">EAN No. 05705243208804</p> <div style="text-align: center;">  <p style="font-size: 8px;">(01)05705243208804(17)170331(10)632574 v.04</p> </div>
<p>...the market unit labelling will include the product code and lot number.</p>	<div style="display: flex; justify-content: space-between; align-items: center; margin-bottom: 10px;"> <div style="font-size: 24px; font-weight: bold;">100</div> <div style="border: 2px solid yellow; padding: 5px; display: inline-block;"> REF 84061181 </div> </div> <p> GB Unometer bag 2000 S Unometerpåse DK Unometerpose 2000 FIN Unometerpussi NO Unometerpose 2000 </p> <div style="display: flex; justify-content: space-between; align-items: center; margin-top: 10px;"> <div style="border: 1px solid black; padding: 2px;"> EC REP UNOMEDICAL A/S Aaholmvej 1-3, Østød 4320 LEJRE DENMARK </div> <div style="text-align: center;">  </div> <div style="font-size: 8px;"> Conva Tec Limited First Avenue Deeside Industrial Park Deeside, Flintshire, CH5 2NU, UK </div> </div> <div style="display: flex; justify-content: space-between; align-items: center; margin-top: 10px;"> <div style="text-align: center;"> <p style="font-size: 8px;">EAN No. 05705243208798</p>  <p style="font-size: 8px;">(01)05705243208798(17)240731(10)632574</p> </div> <div style="text-align: center;">  2019-08  2024-07 </div> <div style="text-align: center;">  </div> </div> <p style="text-align: right; font-size: 8px;">v.05</p>

Attachment 2: Product List

PURPOSE

This attachment is a listing of EU products included in the scope of this Field Safety Notice, which has been organized by Product Code. Within each product code, the lot numbers are organized numerically, smallest to largest.

Please contact your ConvaTec Customer Service Team Member for assistance.

Column Header	Details	What Level of Packaging includes this information?
Product Code	The Product Code is also known as the Catalog Number, or the SKU.	This number is included on the Shipper Label and the Market Unit Label.
Product Description	<p>The Product Description is a description of the finished product from the internal ConvaTec system, SAP. It is a brief, coded description of the product.</p> <p>Utilizing the first line item example, Product Code 84062182, Batch: 156697, UNOMETER BAG 3L BO (50/100) INT:</p> <ul style="list-style-type: none"> • REF: Product code • Batch: Lot number • UNOMETER BAG: UnoMeter™ 500 • BO(50/100): with Bottom Outlet (pack of 50/100) 	A description of the material is on the Shipper Label and on the Market Unit Labeling; however, the description on these labels may not exactly match (verbatim) the Material Description as listed in this attachment.
Lot Number	A traceable number which is specific to the product.	This number is included on the Shipper Label and the Market Unit Label.

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	167273
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	169311
40233182	UNO HANDYVAC CH10 MTL TROC WOBO(1/20)AUS	170793
40233182	UNO HANDYVAC CH10 MTL TROC WOBO(1/20)AUS	172068
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	172069
40097181	HANDYVAC M/700 ML. POSE	172070
40105181	HANDYVAC CH 14 (20/ltk)	172071
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	172072
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	172072
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	172073
40123181	HandyVac m/draen og trocar CH 10	172074
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	172075
40125181	HandyVac m/draen og trocar CH 14	172075
40127181	HandyVac m/draen og trocar CH 18	172076
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	172076
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	172077
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	172078
40101182	UNO HANDYVAC CH10 BAG W/TAP(1/20) INT	172489
22020184	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	173172
22021184	UNO GUEDEL AIRWAY SIZE 3 (50/300) INT	173173
40233182	UNO HANDYVAC CH10 MTL TROC WOBO(1/20)AUS	174559
40101181	HANDYVAC CH 10 (20/ltk)	174560
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	174561
40233182	UNO HANDYVAC CH10 MTL TROC WOBO(1/20)AUS	175474
40233182	UNO HANDYVAC CH10 MTL TROC WOBO 1X20 INT	175474
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	175475
84061182	UNOMETER BAG 2L (100/200) IT	176275
84077182	UNOMETER BAG 2L BO (75/150) IT	176276
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	179215
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	179216
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	179217
40233182	UNO HANDYVAC CH10 MTL TROC WOBO 1X20 INT	179218
40235182	UNO HANDYVAC CH14 MTL TROC WOBO(1/20)AUS	179219
40235182	UNO HANDYVAC CH14 MTL TROC WOBO 1X20 INT	179219
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	179220
40160181	HANDYVAC M/BUNDAFTAPNING CH12	179221
40125181	HandyVac m/draen og trocar CH 14	179222
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	179222

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
84061181	UNOMETER BAG 2L (100/200) INT	179232
22022181	NIELUPUTKI 4 st (4x50 kpl/ltk)	179233
22022181	UNO GUEDEL AIRWAY SIZE 4 (50/200) INT	179233
84010182	UNO CONN CH04-CH18(100/800) INT	179577
84061182	UNOMETER BAG 2L (100/200) IT	179583
84077182	UNOMETER BAG 2L BO (75/150) IT	179591
40233182	UNO HANDYVAC CH10 MTL TROC WOBO 1X20 INT	179924
40235182	UNO HANDYVAC CH14 MTL TROC WOBO 1X20 INT	179925
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	179926
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	179927
40123181	HandyVac m/draen og trocar CH 10	179928
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	179928
40234182	UNO HANDYVAC CH12 MTL TRCR WOBO(1/20)INT	179929
40160181	HANDYVAC M/BUNDAFTAPNING CH12	179930
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	179931
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	179932
40105181	HANDYVAC CH 14 (20/ltk)	179933
22032022	UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA	179949
22034022	UNO GUEDEL AIRWAY SIZE 4 (50/200) ITA	179950
22022181	NIELUPUTKI 4 st (4x50 kpl/ltk)	179955
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	179960
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	181110
40160181	HANDYVAC M/BUNDAFTAPNING CH12	181129
40233182	UNO HANDYVAC CH10 MTL TROC WOBO(1/20)AUS	181132
40105181	HANDYVAC CH 14 (20/ltk)	181140
40125181	HandyVac m/draen og trocar CH 14	181141
22019182	UNO GUEDEL AIRWAY SIZE 1 (100/300) INT	181150
22021182	UNO GUEDEL AIRWAY SIZE 3 (50/300) INT	181151
22023182	UNO GUEDEL AIRWAY SIZE 0 (100/300) INT	181152
22024182	UNO GUEDEL AIRWAY SIZE 00 (100/300) INT	181153
84061181	UNOMETER BAG 2L (100/200) INT	181170
84010181	KATETERANSATS, STERIL	181173
22033022	UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA	181176
40105182	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	181330
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	181895
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	181925
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	181926
40101182	UNO HANDYVAC CH10 BAG W/TAP(1/20) INT	182121
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	182770

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
40233182	UNO HANDYVAC CH10 MTL TROC WOBO 1X20 INT	182774
40235182	UNO HANDYVAC CH14 MTL TROC WOBO(1/20)AUS	182775
40097181	HANDYVAC M/700 ML. POSE	182782
40101182	UNO HANDYVAC CH10 BAG W/TAP(1/20) INT	183061
40227182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS	183089
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	183424
40223182	UNO HANDYVAC CH14 MTL TROC WOBO(1/20)INT	183425
40101181	HANDYVAC M/DRAEN, CH 10/700 ML POSE	183578
40101181	UNO HANDYVAC CH10 BAG W/TAP (1/20) IT	183578
40097181	HANDYVAC M/700 ML. POSE	183595
40105181	HANDYVAC CH 14 (20/ltk)	183596
40123181	HandyVac m/draen og trocar CH 10	183597
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	183598
40125181	HandyVac m/draen og trocar CH 14	183598
40160181	HANDYVAC M/BUNDAFTAPNING CH12	183606
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	183609
22054023	UNO GUEDEL AIRWAY SIZE 5 (1X50PK) IT	183617
40101182	UNO HANDYVAC CH10 BAG W/TAP(1/20) INT	183836
40103182	UNO HANDYVAC CH12 BAG W/TAP(1/20) INT	183837
84010182	UNO CONN CH04-CH18(100/800) INT	183865
84061182	UNOMETER BAG 2L (100/200) IT	183866
84077182	UNOMETER BAG 2L BO (75/150) IT	183871
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROCAR CH14	184213
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	184214
40097181	HANDYVAC M/700 ML. POSE	184230
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	184231
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	184240
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	184241
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	184256
40123181	HandyVac m/draen og trocar CH 10	184257
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	184257
40125181	HandyVac m/draen og trocar CH 14	184258
40127181	HandyVac m/draen og trocar CH 18	184259
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	184259
40160181	HANDYVAC M/BUNDAFTAPNING CH12	184263
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	184273
84061181	UNOMETER BAG 2L (100/200) INT	184292
84010181	KATETERANSATS, STERIL	184300
22054181	TUNGEHOLDER, STR. 5, LYSEBLAA, STERIL	184301

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
22054181	NIELUPUTKI 5 st (50 kpl/ltk)	184301
22020184	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	184770
22021184	UNO GUEDEL AIRWAY SIZE 3 (50/300) INT	184771
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	185141
40234182	UNO HANDYVAC CH12 MTL TRCR WOBO(1/20)INT	185142
40105181	HANDYVAC CH 14 (20/ltk)	185870
84061182	UNOMETER BAG 2L (100/200) IT	186406
84077182	UNOMETER BAG 2L BO (75/150) IT	186412
84062182	UNOMETER BAG 3L BO (50/100) INT	186413
84010182	UNO CONN CH04-CH18(100/800) INT	187037
84062182	UNOMETER BAG 3L BO (50/100) INT	187039
40125181	HandyVac m/draen og trocar CH 14	187521
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	187521
40127181	HandyVac m/draen og trocar CH 18	187522
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	187522
40223183	HANDYVAC + trocarCH14 + BOLSA cerrada 20ud	187531
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	187540
22019182	UNO GUEDEL AIRWAY SIZE 1 (100/300) INT	187567
22021182	UNO GUEDEL AIRWAY SIZE 3 (50/300) INT	187569
22023182	UNO GUEDEL AIRWAY SIZE 0 (100/300) INT	187570
22024182	UNO GUEDEL AIRWAY SIZE 00 (100/300) INT	187571
22033022	UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA	187572
84061181	UNOMETER BAG 2L (100/200) INT	187578
84010181	KATETERANSATS, STERIL	187580
84077181	UNOMETER BAG 3L BO (75/150) INT	187589
84062182	UNOMETER BAG 3L BO (50/100) INT	187980
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	188289
40105181	HANDYVAC CH 14 (20/ltk)	188302
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	188303
22020184	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	188593
84062182	UNOMETER BAG 3L BO (50/100) INT	188599
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	188917
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	188925
22020182	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	188929
22021182	UNO GUEDEL AIRWAY SIZE 3 (50/300) INT	188930
84010182	UNO CONN CH04-CH18(100/800) INT	189330
84061182	UNOMETER BAG 2L (100/200) IT	189331
84077182	UNOMETER BAG 2L BO (75/150) IT	189332
84077182	UNO UNOMETER BAG 2L BO (75/150) IT	189332

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
40125181	HandyVac m/draen og trocar CH 14	189879
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	189879
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	189883
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	190075
84010182	UNO CONN CH04-CH18(100/800) INT	190417
84077182	UNOMETER BAG 2L BO (75/150) IT	190421
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	190815
40223182	UNO HANDYVAC CH14 MTL TROC WOBO(1/20)INT	190822
40123181	HandyVac m/draen og trocar CH 10	190830
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	190830
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	190831
40127181	HandyVac m/draen og trocar CH 18	190832
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	190832
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	190840
40101182	UNO HANDYVAC CH10 BAG W/TAP(1/20) INT	191331
84077182	UNOMETER BAG 2L BO (75/150) IT	191361
84062182	UNOMETER BAG 3L BO (50/100) INT	191363
84077183	UNOMETER BAG 2L BO (75/150) INT	191368
40105181	HANDYVAC CH 14 (20/ltk)	191760
22020184	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	192121
22021184	UNO GUEDEL AIRWAY SIZE 3 (50/300) INT	192122
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	192502
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	192529
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	192529
40105181	HANDYVAC CH 14 (20/ltk)	192529
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	192530
84061181	UNOMETER BAG 2L (100/200) INT	192565
84061182	UNOMETER BAG 2L (100/200) IT	192853
84077182	UNOMETER BAG 2L BO (75/150) IT	192856
84077182	UNO UNOMETER BAG 2L BO (75/150) IT	192856
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	193310
40227182	HANDYVAC U/BUNDAFTAPNING CH12	193311
40227182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS	193311
40123181	HandyVac m/draen og trocar CH 10	193337
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	193337
40105181	HANDYVAC CH 14 (20/ltk)	193370
84010181	KATETERANSATS, STERIL	193381
84010182	UNO CONN CH04-CH18(100/800) INT	193833
22020184	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	193837

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
22021184	UNO GUEDEL AIRWAY SIZE 3 (50/300) INT	193838
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	194258
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	194265
40224182	UNO HANDYVAC CH18 MTL TRCR WOBO(1/20)AUS	194272
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	194275
22033022	UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA	194294
22034022	UNO GUEDEL AIRWAY SIZE 4 (50/200) ITA	194295
84077181	UNOMETER BAG 3L BO (75/150) INT	194302
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	194307
40097181	HANDYVAC M/700 ML. POSE	195031
84077181	UNOMETER BAG 3L BO (75/150) INT	195311
40107182	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	195452
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	195659
22024181	NIELUPUTKI 00 st (100 kpl/ltk)	195667
22024181	UNO GUEDEL AIRWAYS SZ 00 (100/300) INT	195667
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	195671
40101182	UNO HANDYVAC CH10 BAG W/TAP(1/20) INT	196102
40103182	UNO HANDYVAC CH12 BAG W/TAP(1/20) INT	196103
84061182	UNOMETER BAG 2L (100/200) IT	196131
84077182	UNOMETER BAG 2L BO (75/150) IT	196135
84077182	UNO UNOMETER BAG 2L BO (75/150) IT	196135
84077183	UNOMETER BAG 2L BO (75/150) INT	196136
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	196538
22054181	NIELUPUTKI 5 st (50 kpl/ltk)	196546
22054181	UNO GUEDEL AIRWAY SIZE 5 (50/200) INT	196546
84010182	UNO CONN CH04-CH18(100/800) INT	196869
22019184	UNO GUEDEL AIRWAY SIZE 1 (100/300) INT	196871
22020184	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	196872
22021184	UNO GUEDEL AIRWAY SIZE 3 (50/300) INT	196873
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	197370
40097181	HANDYVAC M/700 ML. POSE	197385
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	197386
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	197391
40107183	UNO HANDYVAC CH18 BAG W/TAP (1X1PK) IT	197400
22032022	UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA	197428
22033022	UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA	197429
22022181	NIELUPUTKI 4 st (4x50 kpl/ltk)	197431
84010181	KATETERANSATS, STERIL	197432
22054181	NIELUPUTKI 5 st (50 kpl/ltk)	197435

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
22054181	UNO GUEDEL AIRWAY SIZE 5 (50/200) INT	197435
84077181	UNOMETER BAG 3L BO (75/150) INT	197436
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	197437
22037182	UNO GUEDEL AIRWAY SIZE 000 (100/300) INT	197444
84010182	UNO CONN CH04-CH18(100/800) INT	197686
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	198011
84010181	KATETERANSATS, STERIL	198044
40097181	HANDYVAC M/700 ML. POSE	198991
40127181	HandyVac m/draen og trocar CH 18	198992
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	198992
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	199003
22032022	UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA	199018
22033022	UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA	199019
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin vál	199027
84061182	UNOMETER BAG 2L (100/200) IT	199365
84062182	UNOMETER BAG 3L BO (50/100) INT	199368
84077182	UNOMETER BAG 2L BO (75/150) IT	199369
84077181	UNOMETER BAG 3L BO (75/150) INT	199371
22033022	UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA	199572
40101182	UNO HANDYVAC CH10 BAG W/TAP(1/20) INT	199726
84062182	UNOMETER BAG 3L BO (50/100) INT	199824
22023184	UNO GUEDEL AIRWAY SIZE 0 (100/300) INT	200469
84062182	UNOMETER BAG 3L BO (50/100) INT	200470
40097181	HANDYVAC M/700 ML. POSE	201050
40103181	HANDYVAC M/DRAEN, CH 12/700 ML POSE	201052
40103181	HANDYVAC CH 12 (20/ltk)	201052
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	201053
40123181	HandyVac m/draen og trocar CH 10	201054
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	201054
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	201055
40127181	HandyVac m/draen og trocar CH 18	201055
40160181	HANDYVAC M/BUNDAFTAPNING CH12	201075
40160181	HANDYVAC CH12, TROAKAARI (20 KPL)	201075
40160181	UNO HANDYVAC CH12 MTL TROC BO(1/20) INT	201075
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	201093
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	201113
40105181	HANDYVAC CH 14 (20/ltk)	201113
40223183	HANDYVAC + trocar CH14 + BOLSA cerrada 20ud	201119
22032022	UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA	201126

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
22033022	UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA	201127
84010182	UNO CONN CH04-CH18(100/800) INT	201129
84061181	UNOMETER BAG 2L (100/200) INT	201133
22022181	NIELUPUTKI 4 st (4x50 kpl/ltk)	201134
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin vál	201148
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	201149
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROCER CH14	201955
40123181	HandyVac m/draen og trocar CH 10	201989
40127181	HandyVac m/draen og trocar CH 18	201990
40160181	HANDYVAC M/BUNDAFTAPNING CH12	202031
22019182	UNO GUEDEL AIRWAY SIZE 1 (100/300) INT	202050
84061181	UNOMETER BAG 2L (100/200) INT	202055
84010181	KATETERANSATS, STERIL	202061
40107182	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	202306
84061182	UNOMETER BAG 2L (100/200) IT	202466
84077182	UNOMETER BAG 2L BO (75/150) IT	202470
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	202791
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	202792
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	202793
40228183	HANDYVAC + trocCH16 + BOLSA cerrada 20ud	202798
40222183	HANDYVAC + trocCH10 + BOLSA cerrada 20ud	202799
40227183	HANDYVAC + trocCH12 + BOLSA cerrada 20ud	202800
40227183	UNO HANDYVAC CH12 MTL TROC (1/20) INT	202800
40105182	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	203040
84062182	UNOMETER BAG 3L BO (50/100) INT	203056
40222183	HANDYVAC + trocCH10 + BOLSA cerrada 20ud	203491
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	203492
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin vál	203503
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	203504
84077181	UNOMETER BAG 3L BO (75/150) INT	203711
84010182	UNO CONN CH04-CH18(100/800) INT	204020
22019184	UNO GUEDEL AIRWAY SIZE 1 (100/300) INT	204024
22020184	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	204025
22021184	UNO GUEDEL AIRWAY SIZE 3 (50/300) INT	204026
84062182	UNOMETER BAG 3L BO (50/100) INT	204029
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	204320
40222183	HANDYVAC + trocCH10 + BOLSA cerrada 20ud	204335
40233182	UNO HANDYVAC CH10 MTL TROC WOBO 1X20 INT	204342
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	204349

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	204351
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	204353
22054023	UNO GUEDEL AIRWAY SIZE 5 (50/200) IT	204357
22054023	UNO GUEDEL AIRWAY SIZE 5 (1X50PK) IT	204357
22022181	NIELUPUTKI 4 st (4x50 kpl/ltk)	204358
84010181	KATETERANSATS, STERIL	204359
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin vál	204361
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	204362
22021181	NIELUPUTKI 3 st (6x50 kpl/ltk)	204367
40227182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS	205219
40127181	HandyVac m/draen og trocar CH 18	205242
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	205242
22019184	UNO GUEDEL AIRWAY SIZE 1 (100/300) INT	205691
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	206057
22021184	UNO GUEDEL AIRWAY SIZE 3 (50/300) INT	206072
22022184	UNO GUEDEL AIRWAY SIZE 4 (50/200) INT	206073
84061181	UNOMETER BAG 2L (100/200) INT	206074
22019184	UNO GUEDEL AIRWAY SIZE 1 (100/300) INT	206082
22020184	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	206083
84061182	UNOMETER BAG 2L (100/200) IT	206472
84077182	UNOMETER BAG 2L BO (75/150) IT	206476
84062182	UNOMETER BAG 3L BO (50/100) INT	206477
84077181	UNOMETER BAG 3L BO (75/150) INT	206478
84010182	UNO CONN CH04-CH18(100/800) INT	206968
84062182	UNOMETER BAG 3L BO (50/100) INT	206973
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	207363
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	207375
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	207375
40105181	HANDYVAC CH 14 (20/ltk)	207375
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	207376
40107181	HANDYVAC CH 18 (20/ltk)	207376
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	207376
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	207377
22022182	UNO GUEDEL AIRWAY SIZE 4 (50/200) INT	207426
22021184	UNO GUEDEL AIRWAY SIZE 3 (50/300) INT	207429
22022184	UNO GUEDEL AIRWAY SIZE 4 (50/200) INT	207430
22023184	UNO GUEDEL AIRWAY SIZE 0 (100/300) INT	207431
22055181	NIELUPUTKI 6 st (50 kpl/ltk)	207432
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin vál	207433

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	207434
40103182	UNO HANDYVAC CH12 BAG W/TAP(1/20) INT	207674
22020184	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	207685
22024184	UNO GUEDEL AIRWAY SIZE 00 (100/300) INT	207686
84062182	UNOMETER BAG 3L BO (50/100) INT	207691
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	208231
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	208237
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	208249
84061182	UNOMETER BAG 2L (100/200) IT	208604
84077182	UNOMETER BAG 2L BO (75/150) IT	208605
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	208842
40160181	HANDYVAC M/BUNDAFTAPNING CH12	208852
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	208864
40105181	HANDYVAC CH 14 (20/ltk)	208877
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	208878
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	208878
40107181	HANDYVAC CH 18 (20/ltk)	208878
22021181	TUNGEHOLDER, STR. 3, GUL, STERIL	209032
22021181	NIELUPUTKI 3 st (6x50 kpl/ltk)	209032
22021181	UNO GUEDEL AIRWAY SIZE 3 (50/300) INT	209032
22022181	NIELUPUTKI 4 st (4x50 kpl/ltk)	209033
84062182	UNOMETER BAG 3L BO (50/100) INT	209316
22019181	TUNGEHOLDER, STR. 1, HVID, STERIL	210157
22019181	NIELUPUTKI 1 st (100 kpl/ltk)	210157
22019181	UNO GUEDEL AIRWAY SIZE 1 (100/300) INT	210157
22021181	NIELUPUTKI 3 st (6x50 kpl/ltk)	210158
22022181	NIELUPUTKI 4 st (4x50 kpl/ltk)	210159
40097181	HANDYVAC M/700 ML. POSE	210488
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	210489
40105181	HANDYVAC CH 14 (20/ltk)	210489
40107181	HANDYVAC CH 18 (20/ltk)	210490
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	210508
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	210515
22019181	NIELUPUTKI 1 st (100 kpl/ltk)	210525
22021181	NIELUPUTKI 3 st (6x50 kpl/ltk)	210526
84010181	KATETERANSATS, STERIL	210527
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin vál	210529
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin val	210529
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	210530

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
84010182	UNO CONN CH04-CH18(100/800) INT	210950
84021182	UNO URIMETER BAG 2L (75/300) IT	210953
84062182	UNOMETER BAG 3L BO (50/100) INT	210957
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	211094
40107183	UNO HANDYVAC CH18 BAG W/TAP (1X1PK) IT	211095
22024184	UNO GUEDEL AIRWAY SIZE 00 (100/300) INT	211116
84010182	UNO CONN CH04-CH18(100/800) INT	211197
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROCER CH14	211357
40097181	HANDYVAC M/700 ML. POSE	211375
40125181	HandyVac m/draen og trocar CH 14	211376
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	211376
40160181	HANDYVAC M/BUNDAFTAPNING CH12	211386
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	211414
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	211414
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	211415
40107181	HANDYVAC CH 18 (20/ltk)	211415
22031022	UNO GUEDEL AIRWAY SIZE 1 (100/300) ITA	211434
22022181	NIELUPUTKI 4 st (4x50 kpl/ltk)	211439
22022181	UNO GUEDEL AIRWAY SIZE 4 (50/200) INT	211439
22020184	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	211653
22021184	UNO GUEDEL AIRWAY SIZE 3 (50/300) INT	211654
22022184	UNO GUEDEL AIRWAY SIZE 4 (50/200) INT	211655
22024184	UNO GUEDEL AIRWAY SIZE 00 (100/300) INT	211656
84062182	UNOMETER BAG 3L BO (50/100) INT	211659
22034022	UNO GUEDEL AIRWAY SIZE 4 (50/200) ITA	211917
22023181	NIELUPUTKI 0 st (100 kpl/ltk)	212079
22023181	UNO GUEDEL AIRWAY SIZE 0 (100/300) INT	212079
84062182	UNOMETER BAG 3L BO (50/100) INT	212409
84077181	UNOMETER BAG 3L BO (75/150) INT	212417
40171183	UNO HANDYVAC CH16 BAG W/TAP (1X1PK) IT	212929
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	212939
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	212939
40105181	HANDYVAC CH 14 (20/ltk)	212939
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	212940
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	212940
40107181	HANDYVAC CH 18 (20/ltk)	212940
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	212944
40227182	UNO HANDYVAC CH12 MTL TROC (1/20) IT	213036
40227182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS	213036

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
40125181	HandyVac m/draen og trocar CH 14	213050
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	213050
84061182	UNOMETER BAG 2L (100/200) IT	213305
84077183	UNOMETER BAG 2L BO (75/150) INT	213307
40097181	HANDYVAC M/700 ML. POSE	213804
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	213805
40160181	HANDYVAC M/BUNDAFTAPNING CH12	213809
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin vál	213854
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	213855
22032022	UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA	213856
22033022	UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA	213857
84062182	UNOMETER BAG 3L BO (50/100) INT	214202
40097181	HANDYVAC M/700 ML. POSE	214593
40123181	HandyVac m/draen og trocar CH 10	214594
40125181	HandyVac m/draen og trocar CH 14	214595
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	214595
40160181	HANDYVAC M/BUNDAFTAPNING CH12	214600
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	214609
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	214619
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	215190
40123181	HandyVac m/draen og trocar CH 10	215190
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	215199
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	215203
40107182	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	215426
84077182	UNOMETER BAG 2L BO (75/150) IT	215435
84077182	UNO UNOMETER BAG 2L BO (75/150) IT	215435
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	216421
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	216428
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	216428
40105181	HANDYVAC CH 14 (20/ltk)	216428
40107181	HANDYVAC CH 18 (20/ltk)	216429
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	216432
22034022	UNO GUEDEL AIRWAY SIZE 4 (50/200) ITA	216452
22024184	UNO GUEDEL AIRWAY SIZE 00 (100/300) INT	216457
84062182	UNOMETER BAG 3L BO (50/100) INT	216864
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	217375
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	217386
40107183	UNO HANDYVAC CH18 BAG W/TAP (1X1PK) IT	217387
40103182	UNO HANDYVAC CH12 BAG W/TAP(1/20) INT	217661

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
22020184	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	217681
40223182	UNO HANDYVAC CH14 MTL TROC WOBO(1/20)INT	217723
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	218119
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	218128
40105181	HANDYVAC CH 14 (20/ltk)	218128
40107181	HANDYVAC CH 18 (20/ltk)	218129
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	218130
40123181	HandyVac m/draen og trocar CH 10	218131
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	218131
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	218132
40127181	HandyVac m/draen og trocar CH 18	218133
84010182	UNO CONN CH04-CH18(100/800) INT	218429
84062182	UNOMETER BAG 3L BO (50/100) INT	218441
84077181	UNOMETER BAG 3L BO (75/150) INT	218482
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	218636
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	218641
40107181	HANDYVAC CH 18 (20/ltk)	218645
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	218645
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	218646
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	218653
40103182	UNO HANDYVAC CH12 BAG W/TAP(1/20) INT	219139
84010182	UNO CONN CH04-CH18(100/800) INT	219167
84061182	UNOMETER BAG 2L (100/200) IT	219168
84077182	UNO UNOMETER BAG 2L BO (75/150) IT	219169
84077181	UNOMETER BAG 3L BO (75/150) INT	219175
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	219651
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	219664
40125181	HandyVac m/draen og trocar CH 14	219665
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	219665
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	219690
40105181	HANDYVAC CH 14 (20/ltk)	219690
40107181	HANDYVAC CH 18 (20/ltk)	219691
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	219695
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	219696
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	219698
22020181	NIELUPUTKI 2 st (6x50 kpl/ltk)	219725
22021181	NIELUPUTKI 3 st (6x50 kpl/ltk)	219731
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin val	219733
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	219734

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
84062182	UNOMETER BAG 3L BO (50/100) INT	219982
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	220356
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	220357
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	220362
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	220374
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	220377
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	220798
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	221105
40097181	HANDYVAC M/700 ML. POSE	221123
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	221124
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	221124
40105181	HANDYVAC CH 14 (20/ltk)	221124
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	221125
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	221125
40125181	HandyVac m/draen og trocar CH 14	221126
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	221126
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	221127
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	221144
84077183	UNOMETER BAG 2L BO (75/150) INT	221502
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	221855
40223182	UNO HANDYVAC CH14 MTL TROC WOBO(1/20)INT	222210
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	222548
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	222549
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	222549
40107181	HANDYVAC CH 18 (20/ltk)	222549
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	222554
40097182	UNO HANDYVAC BAG 700ML BO(1/20) INT	222555
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	222569
22021184	UNO GUEDEL AIRWAY SIZE 3 (50/300) INT	223012
22022184	UNO GUEDEL AIRWAY SIZE 4 (50/200) INT	223013
84077182	UNOMETER BAG 2L BO (75/150) IT	223014
84062182	UNOMETER BAG 3L BO (50/100) INT	223018
40123181	HandyVac m/draen og trocar CH 10	223452
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	223452
40125181	HandyVac m/draen og trocar CH 14	223453
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	223453
40160181	HANDYVAC M/BUNDAFTAPNING CH12	223458
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	223459
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	223460

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	223470
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	223479
40234182	UNO HANDYVAC CH12 MTL TRCR WOBO(1/20)INT	223482
40234182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 INT	223482
22032022	UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA	223490
84061181	UNOMETER BAG 2L (100/200) INT	223493
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin val	223496
40223182	UNO HANDYVAC CH14 MTL TROC WOBO(1/20)INT	224343
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	224343
40123181	HandyVac m/draen og trocar CH 10	224364
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	224364
40127181	HandyVac m/draen og trocar CH 18	224365
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	224365
40160181	HANDYVAC M/BUNDAFTAPNING CH12	224374
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	224400
40105181	HANDYVAC CH 14 (20/ltk)	224400
40125181	HandyVac m/draen og trocar CH 14	224401
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	224401
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	224404
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	224408
84010182	UNO CONN CH04-CH18(100/800) INT	224415
40222183	HANDYVAC + trocCH10 + BOLSA cerrada 20ud	224468
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	225138
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	225144
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	225319
84061181	UNOMETER BAG 2L (100/200) INT	225323
84077182	UNOMETER BAG 2L BO (75/150) IT	225721
84062182	UNOMETER BAG 3L BO (50/100) INT	225722
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	226131
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	226132
40107183	UNO HANDYVAC CH18 BAG W/TAP (1X1PK) IT	226133
84010182	UNO CONN CH04-CH18(100/800) INT	226444
22020184	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	226446
22021184	UNO GUEDEL AIRWAY SIZE 3 (50/300) INT	226447
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	226589
40223182	UNO HANDYVAC CH14 MTL TROC WOBO(1/20)INT	226589
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	226876
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	226878
40097181	HANDYVAC M/700 ML. POSE	226893

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
40127181	HandyVac m/draen og trocar CH 18	226895
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	226895
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	226922
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	226923
40107181	HANDYVAC CH 18 (20/ltk)	226923
40234182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 INT	226928
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	226929
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	226932
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	226938
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	226944
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin val	226947
84062182	UNOMETER BAG 3L BO (50/100) INT	227255
84077183	UNOMETER BAG 2L BO (75/150) INT	227257
40097181	HANDYVAC M/700 ML. POSE	227841
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	227842
40105181	HANDYVAC CH 14 (20/ltk)	227842
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	227843
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	227843
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	227844
40160181	HANDYVAC M/BUNDAFTAPNING CH12	227851
84061181	UNOMETER BAG 2L (100/200) INT	227883
40101182	UNO HANDYVAC CH10 BAG W/TAP(1/20) INT	228299
40105182	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	228300
84010182	UNO CONN CH04-CH18(100/800) INT	228333
84061182	UNOMETER BAG 2L (100/200) IT	228335
84077182	UNOMETER BAG 2L BO (75/150) IT	228336
84062182	UNOMETER BAG 3L BO (50/100) INT	228337
40105181	HANDYVAC CH 14 (20/ltk)	228787
40107181	HANDYVAC CH 18 (20/ltk)	228788
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	228791
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin val	228822
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	228824
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	229734
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	229734
40127181	HandyVac m/draen og trocar CH 18	229735
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	229735
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	229739
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	229741
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	229751

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
84061181	UNOMETER BAG 2L (100/200) INT	229767
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin val	229769
84062182	UNOMETER BAG 3L BO (50/100) INT	230224
40160181	HANDYVAC M/BUNDAFTAPNING CH12	230760
40097181	HANDYVAC M/700 ML. POSE	230776
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	230777
40105181	HANDYVAC CH 14 (20/ltk)	230777
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	230778
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	230778
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	230779
84061181	UNOMETER BAG 2L (100/200) INT	230787
40103182	UNO HANDYVAC CH12 BAG W/TAP(1/20) INT	231155
40222183	HANDYVAC + trocCH10 + BOLSA cerrada 20ud	231633
40097181	HANDYVAC M/700 ML. POSE	232455
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	232457
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	232467
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	232468
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	232468
22019184	UNO GUEDEL AIRWAY SIZE 1 (100/300) INT	232760
22020184	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	232761
22021184	UNO GUEDEL AIRWAY SIZE 3 (50/300) INT	232762
22022184	UNO GUEDEL AIRWAY SIZE 4 (50/200) INT	232763
84062182	UNOMETER BAG 3L BO (50/100) INT	232764
40097181	HANDYVAC M/700 ML. POSE	233067
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	233068
22035022	UNO GUEDEL AIRWAY SIZE 0 (100/300) ITA	233096
84061181	UNOMETER BAG 2L (100/200) INT	233098
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin val	233102
84010182	UNO CONN CH04-CH18(100/800) INT	233717
84061182	UNOMETER BAG 2L (100/200) IT	233721
22019184	UNO GUEDEL AIRWAY SIZE 1 (100/300) INT	233722
22020184	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	233723
22021184	UNO GUEDEL AIRWAY SIZE 3 (50/300) INT	233724
22022184	UNO GUEDEL AIRWAY SIZE 4 (50/200) INT	233725
84077182	UNOMETER BAG 2L BO (75/150) IT	233730
84062182	UNOMETER BAG 3L BO (50/100) INT	233733
22021183	UNO GUEDEL AIRWAY SIZE 3 (1X50PK) IT	233740
84077183	UNOMETER BAG 2L BO (75/150) INT	233745
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	233766

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
40107181	HANDYVAC CH 18 (20/ltk)	233766
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	234147
40107183	UNO HANDYVAC CH18 BAG W/TAP (1X1PK) IT	234155
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	234160
40097181	HANDYVAC M/700 ML. POSE	234909
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	234910
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	234911
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	234911
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	234912
84061181	UNOMETER BAG 2L (100/200) INT	234926
22020182	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	234929
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	235130
40223182	UNO HANDYVAC CH14 MTL TROC WOBO 1X20 AUS	235130
84010182	UNO CONN CH04-CH18(100/800) INT	235326
22019184	UNO GUEDEL AIRWAY SIZE 1 (100/300) INT	235328
22020184	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	235329
22021184	UNO GUEDEL AIRWAY SIZE 3 (50/300) INT	235330
22022184	UNO GUEDEL AIRWAY SIZE 4 (50/200) INT	235331
84062182	UNOMETER BAG 3L BO (50/100) INT	235332
22020183	UNO GUEDEL AIRWAY SIZE 2 (1X50PK) IT	235336
22021183	UNO GUEDEL AIRWAY SIZE 3 (1X50PK) IT	235337
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	235729
40097181	HANDYVAC M/700 ML. POSE	236643
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	236644
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	236644
40105181	HANDYVAC CH 14 (20/ltk)	236644
40107181	HANDYVAC CH 18 (20/ltk)	236645
40125181	HandyVac m/draen og trocar CH 14	236646
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	236646
40160181	HANDYVAC M/BUNDAFTAPNING CH12	236652
40223183	HANDYVAC + trocar CH14 + BOLSA cerrada 20ud	236654
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	236676
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	237467
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	237474
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	237476
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	237477
84077183	UNOMETER BAG 2L BO (75/150) INT	300309
84062182	UNOMETER BAG 3L BO (50/100) INT	300333
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	300339

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
84077183	UNOMETER BAG 2L BO (75/150) INT	300339
84061182	UNOMETER BAG 2L (100/200) IT	300353
84077182	UNOMETER BAG 2L BO (75/150) IT	300355
84062182	UNOMETER BAG 3L BO (50/100) INT	300412
22020184	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	300478
22021184	UNO GUEDEL AIRWAY SIZE 3 (50/300) INT	300479
22022184	UNO GUEDEL AIRWAY SIZE 4 (50/200) INT	300480
84062182	UNOMETER BAG 3L BO (50/100) INT	300542
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	300588
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin val	300590
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	300602
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROCER CH14	300608
40097181	HANDYVAC M/700 ML. POSE	300625
40105181	HANDYVAC CH 14 (20/ltk)	300626
40125181	HandyVac m/draen og trocar CH 14	300627
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	300627
40160181	HANDYVAC M/BUNDAFTAPNING CH12	300634
84010181	KATETERANSATS, STERIL	300642
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	300665
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	300678
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	300678
40107181	HANDYVAC CH 18 (20/ltk)	300678
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	300682
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	300686
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin val	300692
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	300693
40105181	HANDYVAC CH 14 (20/ltk)	300705
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	300705
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	300706
40125181	HandyVac m/draen og trocar CH 14	300708
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	300708
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	300808
40123181	HandyVac m/draen og trocar CH 10	300837
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	300837
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	300839
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	300879
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	300892
40105181	HANDYVAC CH 14 (20/ltk)	300892
40123181	HandyVac m/draen og trocar CH 10	300893

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	300893
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	300900
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	300912
40107181	HANDYVAC CH 18 (20/ltk)	300921
40125181	HandyVac m/draen og trocar CH 14	300922
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin val	300924
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	300926
22020184	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	300931
84062182	UNOMETER BAG 3L BO (50/100) INT	300934
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	301369
40105181	HANDYVAC CH 14 (20/ltk)	301393
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	301393
40107181	HANDYVAC CH 18 (20/ltk)	301394
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	301394
40123181	HandyVac m/draen og trocar CH 10	301395
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	301395
40125181	HandyVac m/draen og trocar CH 14	301396
84077181	UNOMETER BAG 3L BO (75/150) INT	301400
40160181	HANDYVAC M/BUNDAFTAPNING CH12	301401
40222183	HANDYVAC + trocCH10 + BOLSA cerrada 20ud	301402
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	301403
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	301405
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	301407
84062182	UNOMETER BAG 3L BO (50/100) INT	301460
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	301489
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	301490
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	301497
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	301498
22020184	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	301551
84062182	UNOMETER BAG 3L BO (50/100) INT	301571
22021183	UNO GUEDEL AIRWAY SIZE 3 (1X50PK) IT	301594
84077181	UNOMETER BAG 3L BO (75/150) INT	301765
40160181	HANDYVAC M/BUNDAFTAPNING CH12	301766
84010181	KATETERANSATS, STERIL	301784
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	301785
40105181	HANDYVAC CH 14 (20/ltk)	301785
40123181	HandyVac m/draen og trocar CH 10	301786
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	301786
40125181	HandyVac m/draen og trocar CH 14	301787

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	301787
84077182	UNOMETER BAG 2L BO (75/150) IT	301841
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	302103
22032022	UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA	302125
40097181	HANDYVAC M/700 ML. POSE	302139
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	302140
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	302140
40107181	HANDYVAC CH 18 (20/ltk)	302140
40222183	HANDYVAC + trocCH10 + BOLSA cerrada 20ud	302144
84062182	UNOMETER BAG 3L BO (50/100) INT	302234
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	302277
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	302286
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin val	302301
40227183	UNO HANDYVAC CH12 MTL TROC (1/20) INT	302448
40107181	HANDYVAC CH 18 (20/ltk)	302527
40123181	HandyVac m/draen og trocar CH 10	302528
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	302528
40160181	HANDYVAC M/BUNDAFTAPNING CH12	302533
22020182	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	302541
22020184	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	302557
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	302560
40105181	HANDYVAC CH 14 (20/ltk)	302560
84062182	UNOMETER BAG 3L BO (50/100) INT	302587
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	302760
84077183	UNOMETER BAG 2L BO (75/150) INT	302760
40222183	HANDYVAC + trocCH10 + BOLSA cerrada 20ud	302786
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	302787
40107181	HANDYVAC CH 18 (20/ltk)	302797
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	302798
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	302809
40227182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS	302810
84077181	UNOMETER BAG 3L BO (75/150) INT	302855
22022183	UNO GUEDEL AIRWAY SIZE 4 (50/200) INT	302864
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	302945
40097181	HANDYVAC M/700 ML. POSE	302971
40103181	HANDYVAC M/DRAEN, CH 12/700 ML POSE	302972
40103181	HANDYVAC CH 12 (20/ltk)	302972
40123181	HandyVac m/draen og trocar CH 10	302973
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	302973

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
40125181	HandyVac m/draen og trocar CH 14	302974
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	302974
40127181	HandyVac m/draen og trocar CH 18	302975
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	302975
84077181	UNOMETER BAG 3L BO (75/150) INT	302989
40160181	HANDYVAC M/BUNDAFTAPNING CH12	302990
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	303009
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	303009
40105181	HANDYVAC CH 14 (20/ltk)	303009
40234182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 INT	303011
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	303012
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	303014
40101182	UNO HANDYVAC CH10 BAG W/TAP(1/20) INT	303176
40105182	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	303177
84077181	UNOMETER BAG 3L BO (75/150) INT	303178
84061182	UNOMETER BAG 2L (100/200) IT	303297
84077182	UNOMETER BAG 2L BO (75/150) IT	303321
84062182	UNOMETER BAG 3L BO (50/100) INT	303343
84021182	UNO URIMETER BAG 2L (75/300) IT	303382
84010183	UNO CONN CH04-CH18(100/800) IT	303383
40125181	HandyVac m/draen og trocar CH 14	303403
22032022	UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA	303407
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin val	303410
22034022	UNO GUEDEL AIRWAY SIZE 4 (50/200) ITA	303414
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	303427
40227183	HANDYVAC + trocCH12 + BOLSA cerrada 20ud	303457
84010182	UNO CONN CH04-CH18(100/800) INT	303471
84062182	UNOMETER BAG 3L BO (50/100) INT	303501
84010182	UNO CONN CH04-CH18(100/800) INT	303530
40101181	HANDYVAC M/DRAEN, CH 10/700 ML POSE	303585
40101181	UNO HANDYVAC CH10 BAG W/TAP (1/20) IT	303585
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	303586
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	303586
40107181	HANDYVAC CH 18 (20/ltk)	303586
40107182	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	303683
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	303686
40227182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS	303692
22032022	UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA	303761
84010181	KATETERANSATS, STERIL	303769

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
22020182	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	303771
22033022	UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA	303772
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	303781
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	303784
84077182	UNOMETER BAG 2L BO (75/150) IT	303956
84010181	KATETERANSATS, STERIL	304094
40097181	HANDYVAC M/700 ML. POSE	304098
40107181	HANDYVAC CH 18 (20/ltk)	304099
40127181	HandyVac m/draen og trocar CH 18	304100
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	304100
40125181	HandyVac m/draen og trocar CH 14	304123
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	304123
40125181	HANDYVAC CH14, TROAKAARI (20 KPL)	304123
84062182	UNOMETER BAG 3L BO (50/100) INT	304184
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	304271
22055181	NIELUPUTKI 6 st (50 kpl/ltk)	304372
84010182	UNO CONN CH04-CH18(100/800) INT	304391
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	304401
84010181	KATETERANSATS, STERIL	304442
40105181	HANDYVAC CH 14 (20/ltk)	304448
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	304449
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	304449
40107181	HANDYVAC CH 18 (20/ltk)	304449
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	304450
40125181	HandyVac m/draen og trocar CH 14	304451
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	304451
40127181	HandyVac m/draen og trocar CH 18	304452
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	304452
40234182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 INT	304465
40160181	HANDYVAC M/BUNDAFTAPNING CH12	304466
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin val	304468
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	304470
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	304471
84062182	UNOMETER BAG 3L BO (50/100) INT	304623
40103182	UNO HANDYVAC CH12 BAG W/TAP(1/20) INT	304702
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin val	304734
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin val	304737
40222183	HANDYVAC + trocCH10 + BOLSA cerrada 20ud	304738
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	304741

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	304763
84062182	UNOMETER BAG 3L BO (50/100) INT	304798
84077181	UNOMETER BAG 3L BO (75/150) INT	305018
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	305035
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	305035
40107181	HANDYVAC CH 18 (20/ltk)	305035
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	305036
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	305037
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	305053
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	305054
40101181	HANDYVAC M/DRAEN, CH 10/700 ML POSE	305409
40101181	UNO HANDYVAC CH10 BAG W/TAP (1/20) IT	305409
40101181	HANDYVAC CH 10 (20/ltk)	305409
40101182	UNO HANDYVAC CH10 BAG W/TAP(1/20) INT	305410
40103181	HANDYVAC M/DRAEN, CH 12/700 ML POSE	305411
40103181	HANDYVAC CH 12 (20/ltk)	305411
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	305412
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	305412
40105181	HANDYVAC CH 14 (20/ltk)	305412
40107181	HANDYVAC CH 18 (20/ltk)	305413
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	305414
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	305415
40125181	HandyVac m/draen og trocar CH 14	305416
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	305416
40125181	HANDYVAC CH14, TROAKAARI (20 KPL)	305416
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	305417
40127181	HandyVac m/draen og trocar CH 18	305418
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	305418
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	305419
40160181	HANDYVAC M/BUNDAFTAPNING CH12	305425
84010181	KATETERANSATS, STERIL	305438
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin val	305441
84062182	UNOMETER BAG 3L BO (50/100) INT	305444
84077181	UNOMETER BAG 3L BO (75/150) INT	305446
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	305447
40101181	HANDYVAC CH 10 (20/ltk)	305737
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	305738
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	305738
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	305739

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	305739
40107181	HANDYVAC CH 18 (20/ltk)	305739
40125181	HandyVac m/draen og trocar CH 14	305740
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	305740
40127181	HandyVac m/draen og trocar CH 18	305741
40160181	HANDYVAC M/BUNDAFTAPNING CH12	305745
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROCAR CH14	305747
84010181	KATETERANSATS, STERIL	305756
84077181	UNOMETER BAG 3L BO (75/150) INT	305760
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	305766
40223182	UNO HANDYVAC CH14 MTL TROC WOBO 1X20 AUS	306028
84061182	UNOMETER BAG 2L (100/200) IT	306059
84062182	UNOMETER BAG 3L BO (50/100) INT	306060
84077182	UNOMETER BAG 2L BO (75/150) IT	306061
40101182	UNO HANDYVAC CH10 BAG W/TAP(1/20) INT	306194
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	306195
84010181	KATETERANSATS, STERIL	306197
40234182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 INT	306221
40101181	UNO HANDYVAC CH10 BAG W/TAP (1/20) IT	306346
40101181	HANDYVAC CH 10 (20/ltk)	306346
40127181	HandyVac m/draen og trocar CH 18	306347
84010181	KATETERANSATS, STERIL	306358
40101181	HANDYVAC CH 10 (20/ltk)	306452
40101181	UNO HANDYVAC CH10 BAG W/TAP (1/20) IT	306452
40101181	UNO HANDYVAC CH10 BAG W/TAP (1/20) IT	306635
40101181	HANDYVAC CH 10 (20/ltk)	306635
84077182	UNOMETER BAG 2L BO (75/150) IT	306644
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	306645
40101181	UNO HANDYVAC CH10 BAG W/TAP (1/20) IT	306857
40101181	HANDYVAC CH 10 (20/ltk)	306857
40123181	HandyVac m/draen og trocar CH 10	306858
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	306858
40123181	HANDYVAC CH10, TROAKAARI (20 KPL)	306858
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	306859
40097181	HANDYVAC M/700 ML. POSE	306874
84010182	UNO CONN CH04-CH18(100/800) INT	306875
84077181	UNOMETER BAG 3L BO (75/150) INT	306876
40097181	HANDYVAC M/700 ML. POSE	307109
40101181	UNO HANDYVAC CH10 BAG W/TAP (1/20) IT	307110

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
40101181	HANDYVAC CH 10 (20/ltk)	307110
40101182	UNO HANDYVAC CH10 BAG W/TAP(1/20) INT	307111
40103182	UNO HANDYVAC CH12 BAG W/TAP(1/20) INT	307112
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	307113
40105181	HANDYVAC CH 14 (20/ltk)	307113
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	307114
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	307115
40160181	HANDYVAC M/BUNDAFTAPNING CH12	307117
84062182	UNOMETER BAG 3L BO (50/100) INT	307131
84077181	UNOMETER BAG 3L BO (75/150) INT	307133
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	307134
40101181	UNO HANDYVAC CH10 BAG W/TAP (1/20) IT	307353
40101181	HANDYVAC CH 10 (20/ltk)	307353
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	307354
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	307354
40105181	HANDYVAC CH 14 (20/ltk)	307354
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	307355
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	307355
40107181	HANDYVAC CH 18 (20/ltk)	307355
84061182	UNOMETER BAG 2L (100/200) IT	307364
84077182	UNOMETER BAG 2L BO (75/150) IT	307366
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	307367
40105182	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	307379
40160181	HANDYVAC M/BUNDAFTAPNING CH12	307411
40234182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 INT	307412
40097181	HANDYVAC M/700 ML. POSE	307561
40123181	HandyVac m/draen og trocar CH 10	307562
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	307562
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	307563
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	307564
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	307565
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	307568
40235182	UNO HANDYVAC CH14 MTL TROC WOBO 1X20 INT	307569
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	307581
22033022	UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA	307582
40097181	HANDYVAC M/700 ML. POSE	307757
40123181	HandyVac m/draen og trocar CH 10	307758
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	307758
40125181	HandyVac m/draen og trocar CH 14	307759

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	307759
40125181	HANDYVAC CH14, TROAKAARI (20 KPL)	307759
84010181	KATETERANSATS, STERIL	307774
84062182	UNOMETER BAG 3L BO (50/100) INT	307778
40101182	UNO HANDYVAC CH10 BAG W/TAP(1/20) INT	307791
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	307792
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	307794
40222183	HANDYVAC + trocCH10 + BOLSA cerrada 20ud	308155
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	308157
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin val	308169
84077182	UNOMETER BAG 2L BO (75/150) IT	308173
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	308174
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	308181
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	308181
40235182	UNO HANDYVAC CH14 MTL TROC WOBO 1X20 INT	308187
40227182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS	308309
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	308385
40107181	HANDYVAC CH 18 (20/ltk)	308385
84062182	UNOMETER BAG 3L BO (50/100) INT	308392
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	308446
84010182	UNO CONN CH04-CH18(100/800) INT	308449
84077181	UNOMETER BAG 3L BO (75/150) INT	308451
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	308452
84077183	UNOMETER BAG 2L BO (75/150) INT	308452
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	308454
40101181	UNO HANDYVAC CH10 BAG W/TAP (1/20) IT	308461
40101181	HANDYVAC CH 10 (20/ltk)	308461
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	308626
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	308627
40234182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 INT	308628
40227182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS	308664
40097181	HANDYVAC M/700 ML. POSE	308807
40101182	UNO HANDYVAC CH10 BAG W/TAP(1/20) INT	308808
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	308809
40105181	HANDYVAC CH 14 (20/ltk)	308809
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	308810
40127181	HANDYVAC CH18, TROAKAARI (20 KPL)	308811
40160181	HANDYVAC M/BUNDAFTAPNING CH12	308814
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	308816

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
40223182	UNO HANDYVAC CH14 MTL TROC WOBO 1X20 AUS	308816
84010181	KATETERANSATS, STERIL	308822
84077181	UNOMETER BAG 3L BO (75/150) INT	308825
22034022	UNO GUEDEL AIRWAY SIZE 4 (50/200) ITA	308837
40227182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS	308846
40123181	HandyVac m/draen og trocar CH 10	308972
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	308972
84010182	UNO CONN CH04-CH18(100/800) INT	308976
84062182	UNOMETER BAG 3L BO (50/100) INT	308979
40125181	HandyVac m/draen og trocar CH 14	309003
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	309003
22032022	UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA	309129
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	309181
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	309181
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	309182
40107182	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	309183
40223183	HANDYVAC + trocarCH14 + BOLSA cerrada 20ud	309185
40227182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS	309186
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	309193
84077183	UNOMETER BAG 2L BO (75/150) INT	309193
40097181	HANDYVAC M/700 ML. POSE	309226
40101181	HANDYVAC M/DRAEN, CH 10/700 ML POSE	309227
40101181	UNO HANDYVAC CH10 BAG W/TAP (1/20) IT	309227
40101181	HANDYVAC CH 10 (20/ltk)	309227
40103181	HANDYVAC M/DRAEN, CH 12/700 ML POSE	309228
40103181	HANDYVAC CH 12 (20/ltk)	309228
40127181	HandyVac m/draen og trocar CH 18	309229
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	309229
84010181	KATETERANSATS, STERIL	309236
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	309422
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	309423
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	309424
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	309425
84010182	UNO CONN CH04-CH18(100/800) INT	309454
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	309563
45 951 1S-50	UNO PLEURA PUNCTURE (1/20) INT	309567
22032022	UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA	309599
40097181	HANDYVAC M/700 ML. POSE	309665
40101182	UNO HANDYVAC CH10 BAG W/TAP(1/20) INT	309666

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	309667
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	309667
40105181	HANDYVAC CH 14 (20/ltk)	309667
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	309668
40107181	HANDYVAC CH 18 (20/ltk)	309668
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	309669
40125181	HandyVac m/draen og trocar CH 14	309670
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	309670
40127181	HandyVac m/draen og trocar CH 18	309671
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	309671
40127181	HANDYVAC CH18, TROAKAARI (20 KPL)	309671
40160181	HANDYVAC M/BUNDAFTAPNING CH12	309674
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	309676
40223182	UNO HANDYVAC CH14 MTL TROC WOBO 1X20 AUS	309676
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	309677
40235182	UNO HANDYVAC CH14 MTL TROC WOBO 1X20 INT	309680
84010181	KATETERANSATS, STERIL	309691
84062182	UNOMETER BAG 3L BO (50/100) INT	309694
84077181	UNOMETER BAG 3L BO (75/150) INT	309696
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	309697
84077183	UNOMETER BAG 2L BO (75/150) INT	309697
40097181	HANDYVAC M/700 ML. POSE	309885
40125181	HandyVac m/draen og trocar CH 14	309886
40127181	HandyVac m/draen og trocar CH 18	309887
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	309887
40127181	HANDYVAC CH18, TROAKAARI (20 KPL)	309887
40101182	UNO HANDYVAC CH10 BAG W/TAP(1/20) INT	309916
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	309917
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	309918
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	309919
84010181	KATETERANSATS, STERIL	309923
22033022	UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA	309932
22034022	UNO GUEDEL AIRWAY SIZE 4 (50/200) ITA	309933
40097181	HANDYVAC M/700 ML. POSE	310091
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	310092
40105181	HANDYVAC CH 14 (20/ltk)	310092
40107181	HANDYVAC CH 18 (20/ltk)	310093
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	310094
40125181	HandyVac m/draen og trocar CH 14	310095

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	310096
40222183	HANDYVAC + trocCH10 + BOLSA cerrada 20ud	310099
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	310100
84010182	UNO CONN CH04-CH18(100/800) INT	310112
84062182	UNOMETER BAG 3L BO (50/100) INT	310117
84077181	UNOMETER BAG 3L BO (75/150) INT	310119
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	310120
40184182	UNO HANDYVAC CH16 MTL TROCAR BO(1/20) IT	310144
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	310391
40125181	HandyVac m/draen og trocar CH 14	310392
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	310392
40125181	HANDYVAC CH14, TROAKAARI (20 KPL)	310392
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	310403
40097181	HANDYVAC M/700 ML. POSE	310405
40101182	UNO HANDYVAC CH10 BAG W/TAP(1/20) INT	310406
40127181	HandyVac m/draen og trocar CH 18	310423
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	310423
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	310578
40103182	UNO HANDYVAC CH12 BAG W/TAP(1/20) INT	310597
40123181	HandyVac m/draen og trocar CH 10	310598
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	310598
84062182	UNOMETER BAG 3L BO (50/100) INT	310603
22020184	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	310605
22021184	UNO GUEDEL AIRWAY SIZE 3 (50/300) INT	310606
40097181	HANDYVAC M/700 ML. POSE	310786
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	310787
40105181	HANDYVAC CH 14 (20/ltk)	310787
40123181	HandyVac m/draen og trocar CH 10	310788
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	310788
40125181	HandyVac m/draen og trocar CH 14	310789
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	310789
40125181	HANDYVAC CH14, TROAKAARI (20 KPL)	310789
40127181	HandyVac m/draen og trocar CH 18	310790
40127181	HANDYVAC CH18, TROAKAARI (20 KPL)	310790
40160181	HANDYVAC M/BUNDAFTAPNING CH12	310793
84010181	KATETERANSATS, STERIL	310805
84062182	UNOMETER BAG 3L BO (50/100) INT	310824
84077181	UNOMETER BAG 3L BO (75/150) INT	310825
22032022	UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA	310976

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
22033022	UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA	310977
40097181	HANDYVAC M/700 ML. POSE	311012
40103181	HANDYVAC CH 12 (20/ltk)	311013
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	311014
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin val	311019
84061182	UNOMETER BAG 2L (100/200) IT	311022
84077182	UNOMETER BAG 2L BO (75/150) IT	311023
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	311032
40107181	HANDYVAC CH 18 (20/ltk)	311186
84021182	UNO URIMETER BAG 2L (75/300) IT	311190
84062182	UNOMETER BAG 3L BO (50/100) INT	311201
22032022	UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA	311326
22033022	UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA	311327
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	311367
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	311367
40105181	HANDYVAC CH 14 (20/ltk)	311367
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	311368
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	311368
40107181	HANDYVAC CH 18 (20/ltk)	311368
40160181	HANDYVAC M/BUNDAFTAPNING CH12	311369
40101181	HANDYVAC M/DRAEN, CH 10/700 ML POSE	311437
40101181	UNO HANDYVAC CH10 BAG W/TAP (1/20) IT	311437
40101181	HANDYVAC CH 10 (20/ltk)	311437
84077182	UNOMETER BAG 2L BO (75/150) IT	311571
22020184	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	311680
22021184	UNO GUEDEL AIRWAY SIZE 3 (50/300) INT	311681
22022184	UNO GUEDEL AIRWAY SIZE 4 (50/200) INT	311682
40222183	HANDYVAC + trocCH10 + BOLSA cerrada 20ud	311699
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	311700
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin val	311702
40101181	HANDYVAC M/DRAEN, CH 10/700 ML POSE	311717
40101181	UNO HANDYVAC CH10 BAG W/TAP (1/20) IT	311717
40101181	HANDYVAC CH 10 (20/ltk)	311717
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	311718
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	311735
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	311737
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	312059
40105181	HANDYVAC CH 14 (20/ltk)	312059
40107182	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	312060

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	312062
40127181	HandyVac m/draen og trocar CH 18	312063
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	312063
40227182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS	312067
40235182	UNO HANDYVAC CH14 MTL TROC WOBO 1X20 INT	312069
40103181	HANDYVAC CH 12 (20/ltk)	312108
84010182	UNO CONN CH04-CH18(100/800) INT	312205
40125181	HandyVac m/draen og trocar CH 14	312251
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	312251
40125181	HANDYVAC CH14, TROAKAARI (20 KPL)	312251
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin val	312360
22021183	UNO GUEDEL AIRWAY SIZE 3 (1X50PK) IT	312370
40101182	UNO HANDYVAC CH10 BAG W/TAP(1/20) INT	312386
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	312390
40171183	UNO HANDYVAC CH16 BAG W/TAP (1X1PK) IT	312475
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	312498
40223182	UNO HANDYVAC CH14 MTL TROC WOBO 1X20 AUS	312498
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin val	312501
22020183	UNO GUEDEL AIRWAY SIZE 2 (1X50PK) IT	312629
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	312715
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	312715
40105181	HANDYVAC CH 14 (20/ltk)	312715
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	312716
40107181	HANDYVAC CH 18 (20/ltk)	312716
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	312717
84062182	UNOMETER BAG 3L BO (50/100) INT	312733
22033022	UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA	312870
84010182	UNO CONN CH04-CH18(100/800) INT	312884
84010181	KATETERANSATS, STERIL	312891
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	313032
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	313032
40107181	HANDYVAC CH 18 (20/ltk)	313032
84010181	KATETERANSATS, STERIL	313038
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	313041
22021183	UNO GUEDEL AIRWAY SIZE 3 (1X50PK) IT	313054
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	313311
40107181	HANDYVAC CH 18 (20/ltk)	313343
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	313470
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	313470

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
84010181	KATETERANSATS, STERIL	313479
84061182	UNOMETER BAG 2L (100/200) IT	313481
84077182	UNOMETER BAG 2L BO (75/150) IT	313482
84077182	UNO UNOMETER BAG 2L BO (75/150) IT	313482
40125181	HandyVac m/draen og trocar CH 14	313500
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	313500
84010182	UNO CONN CH04-CH18(100/800) INT	313798
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	313799
40097181	HANDYVAC M/700 ML. POSE	313825
84077182	UNOMETER BAG 2L BO (75/150) IT	313828
84077182	UNO UNOMETER BAG 2L BO (75/150) IT	313828
40101181	UNO HANDYVAC CH10 BAG W/TAP (1/20) IT	313845
40101181	HANDYVAC CH 10 (20/ltk)	313845
40101181	HANDYVAC M/DRAEN, CH 10/700 ML POSE	313845
40103181	HANDYVAC M/DRAEN, CH 12/700 ML POSE	313846
40103181	HANDYVAC CH 12 (20/ltk)	313846
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	313847
40125181	HandyVac m/draen og trocar CH 14	313848
40127181	HandyVac m/draen og trocar CH 18	313849
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	313849
40127181	HANDYVAC CH18, TROAKAARI (20 KPL)	313849
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	314075
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	314075
40107181	HANDYVAC CH 18 (20/ltk)	314075
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	314077
40103182	UNO HANDYVAC CH12 BAG W/TAP(1/20) INT	314099
84062182	UNOMETER BAG 3L BO (50/100) INT	314100
22020182	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	314111
22020183	UNO GUEDEL AIRWAY SIZE 2 (1X50PK) IT	314112
22054182	UNO GUEDEL AIRWAY SIZE 5 (50/200) INT	314114
40097181	HANDYVAC M/700 ML. POSE	314271
40123181	HandyVac m/draen og trocar CH 10	314274
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	314274
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	314290
84077183	UNOMETER BAG 2L BO (75/150) INT	314290
40101182	UNO HANDYVAC CH10 BAG W/TAP(1/20) INT	314459
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	314583
22032022	UNO GUEDEL AIRWAY SIZE 2 (50/300) ITA	314764
40099181	HANDYVAC CH 8 (20/ltk)	314863

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
84010183	UNO CONN CH04-CH18(100/800) IT	314883
40105182	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	315047
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	315170
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	315170
84077181	UNOMETER BAG 3L BO (75/150) INT	315188
40227182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS	315244
40107182	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	315397
40125181	HandyVac m/draen og trocar CH 14	315541
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	315541
84010181	KATETERANSATS, STERIL	315554
84062182	UNOMETER BAG 3L BO (50/100) INT	315597
84021183	UNOMETER 500 BOLSAS 1'5L. 75uds (sin val	315672
84021183	UNO URIMETER BAG 2L (75/300) INT	315672
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	315682
84077183	UNOMETER BAG 2L BO (75/150) INT	315682
22020183	UNO GUEDEL AIRWAY SIZE 2 (1X50PK) IT	315745
22021183	UNO GUEDEL AIRWAY SIZE 3 (1X50PK) IT	315746
84021182	UNO URIMETER BAG 2L (75/300) IT	218437
84061182	UNOMETER BAG 2L (100/200) IT	316293
84021182	UNO URIMETER BAG 2L (75/300) IT	304154
84021182	UNO URIMETER BAG 2L (75/300) IT	308337
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	316306
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	316306
40105181	HANDYVAC CH 14 (20/ltk)	316306
40160181	HANDYVAC M/BUNDAFTAPNING CH12	316308
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	316442
22033022	UNO GUEDEL AIRWAY SIZE 3 (50/300) ITA	316550
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	316569
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	316569
40105181	HANDYVAC CH 14 (20/ltk)	316569
84062182	UNOMETER BAG 3L BO (50/100) INT	316577
40125181	HandyVac m/draen og trocar CH 14	316769
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	316769
40125181	HANDYVAC CH14, TROAKAARI (20 KPL)	316769
40222183	HANDYVAC + trocCH10 + BOLSA cerrada 20ud	316772
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	316773
40234182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 INT	316775
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	316817
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	316818

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
40160181	HANDYVAC M/BUNDAFTAPNING CH12	316820
40160181	UNO HANDYVAC CH12 MTL TROC BO(1/20) INT	316820
84077182	UNOMETER BAG 2L BO (75/150) IT	316825
84077182	UNO UNOMETER BAG 2L BO (75/150) IT	316825
40097181	HANDYVAC M/700 ML. POSE	316986
84010182	UNO CONN CH04-CH18(100/800) INT	317016
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	317117
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	317117
40107181	HANDYVAC CH 18 (20/ltk)	317117
40125181	HandyVac m/draen og trocar CH 14	317118
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	317118
40125181	HANDYVAC CH14, TROAKAARI (20 KPL)	317118
40097181	HANDYVAC M/700 ML. POSE	317308
40223183	HANDYVAC + trocarCH14 + BOLSA cerrada 20ud	317504
84010181	KATETERANSATS, STERIL	317508
22020183	UNO GUEDEL AIRWAY SIZE 2 (50/300) IT	317716
22020183	UNO GUEDEL AIRWAY SIZE 2 (1X50PK) IT	317716
45 951 1S-50	UNO PLEURA PUNCTURE (1/20) INT	317791
84062182	UNOMETER BAG 3L BO (50/100) INT	317991
84077182	UNOMETER BAG 2L BO (75/150) IT	317992
84077182	UNO UNOMETER BAG 2L BO (75/150) IT	317992
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	318069
40099181	HANDYVAC CH 8 (20/ltk)	318168
40101182	UNO HANDYVAC CH10 BAG W/TAP(1/20) INT	318456
40127181	HandyVac m/draen og trocar CH 18	318525
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	318525
22022183	UNO GUEDEL AIRWAY SIZE 4 (50/200) INT	318604
22021182	UNO GUEDEL AIRWAY SIZE 3 (50/300) INT	318758
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	318781
84061182	UNOMETER BAG 2L (100/200) IT	319365
22022183	UNO GUEDEL AIRWAY SIZE 4 (50/200) INT	319385
84062182	UNOMETER BAG 3L BO (50/100) INT	319449
84010181	KATETERANSATS, STERIL	319796
40103182	UNO HANDYVAC CH12 BAG W/TAP(1/20) INT	319816
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	319817
40123181	HandyVac m/draen og trocar CH 10	319818
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	319818
40123181	HANDYVAC CH10, TROAKAARI (20 KPL)	319818
40123181	HandyVac m/draen og trocar CH 10	319819

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
40123181	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	319819
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	319820
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	319821
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	319822
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	319825
40223182	UNO HANDYVAC CH14 MTL TROC WOBO 1X20 AUS	319825
40223183	HANDYVAC + trocCH14 + BOLSA cerrada 20ud	319826
40227182	HANDYVAC U/BUNDAFTAPNING CH12	319828
40227182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS	319828
40235182	UNO HANDYVAC CH14 MTL TROC WOBO 1X20 INT	319829
40101181	HANDYVAC M/DRAEN, CH 10/700 ML POSE	319830
40101181	UNO HANDYVAC CH10 BAG W/TAP (1/20) IT	319830
40101181	HANDYVAC CH 10 (20/ltk)	319830
40097181	HANDYVAC M/700 ML. POSE	319962
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	320024
40222183	HANDYVAC + trocCH10 + BOLSA cerrada 20ud	320284
84077182	UNO UNOMETER BAG 2L BO (75/150) IT	320287
22020183	UNO GUEDEL AIRWAY SIZE 2 (1X50PK) IT	320337
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	320347
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	320437
84062182	UNOMETER BAG 3L BO (50/100) INT	320558
40103182	UNO HANDYVAC CH12 BAG W/TAP(1/20) INT	320572
40101182	UNO HANDYVAC CH10 BAG W/TAP(1/20) INT	320582
40127181	HandyVac m/draen og trocar CH 18	320679
40127181	HANDYVAC CH18, TROAKAARI (20 KPL)	320679
40103181	HANDYVAC M/DRAEN, CH 12/700 ML POSE	320728
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	320835
84010182	UNO CONN CH04-CH18(100/800) INT	320846
22037182	UNO GUEDEL AIRWAY SIZE 000 (100/300) INT	320897
40105182	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	320994
84061182	UNOMETER BAG 2L (100/200) IT	320995
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	321027
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	321027
40105181	HANDYVAC CH 14 (20/ltk)	321027
40107182	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	321028
84062182	UNOMETER BAG 3L BO (50/100) INT	321163
40107181	HANDYVAC CH 18 (20/ltk)	321211
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	321212
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	321213

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	321214
40127181	HandyVac m/draen og trocar CH 18	321214
40127181	HANDYVAC CH18, TROAKAARI (20 KPL)	321214
40127182	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	321215
22020182	UNO GUEDEL AIRWAY SIZE 2 (50/300) INT	321407
22021182	UNO GUEDEL AIRWAY SIZE 3 (50/300) INT	321408
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	321487
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	321487
40105181	HANDYVAC CH 14 (20/ltk)	321487
40127181	HandyVac m/draen og trocar CH 18	321488
40127181	HANDYVAC CH18, TROAKAARI (20 KPL)	321488
84061182	UNOMETER BAG 2L (100/200) IT	321505
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	321528
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	321528
40107181	HANDYVAC CH 18 (20/ltk)	321528
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	321640
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	321641
40235182	UNO HANDYVAC CH14 MTL TROC WOBO 1X20 INT	321644
40097181	HANDYVAC M/700 ML. POSE	321690
22020183	UNO GUEDEL AIRWAY SIZE 2 (50/300) IT	321694
84010182	UNO CONN CH04-CH18(100/800) INT	321887
84077182	UNOMETER BAG 2L BO (75/150) IT	321924
84077182	UNO UNOMETER BAG 2L BO (75/150) IT	321924
84077181	UNOMETER BAG 3L BO (75/150) INT	322228
40160181	HANDYVAC M/BUNDAFTAPNING CH12	322411
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	322418
40127181	HandyVac m/draen og trocar CH 18	322428
40127181	UNO HANDYVAC CH18 MTL TROC BO(1/20) INT	322428
40127181	HANDYVAC CH18, TROAKAARI (20 KPL)	322428
40160181	HANDYVAC M/BUNDAFTAPNING CH12	322886
40227182	HANDYVAC U/BUNDAFTAPNING CH12	322891
40227182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS	322891
84077182	UNO UNOMETER BAG 2L BO (75/150) IT	322898
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	322917
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROC CAR CH14	322941
40223182	UNO HANDYVAC CH14 MTL TROC WOBO 1X20 AUS	322941
84062182	UNOMETER BAG 3L BO (50/100) INT	322972
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	323025
40097181	HANDYVAC M/700 ML. POSE	323288

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
84061182	UNOMETER BAG 2L (100/200) IT	323383
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	323552
84077182	UNOMETER BAG 2L BO (75/150) IT	323599
84077182	UNO UNOMETER BAG 2L BO (75/150) IT	323599
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	323702
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	323702
40105181	HANDYVAC CH 14 (20/ltk)	323702
84077181	UNOMETER BAG 3L BO (75/150) INT	323705
40227183	HANDYVAC + trocCH12 + BOLSA cerrada 20ud	323710
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROCAR CH14	323893
40223182	UNO HANDYVAC CH14 MTL TROC WOBO 1X20 AUS	323893
84062182	UNOMETER BAG 3L BO (50/100) INT	323897
40227182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS	324189
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	324192
40107182	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	324205
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	324206
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	324363
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	324363
40105181	HANDYVAC CH 14 (20/ltk)	324363
40125182	UNO HANDYVAC CH14 MTL TROC BO(1/20) EUR	324373
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	324407
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	324459
40125181	HandyVac m/draen og trocar CH 14	324619
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	324619
40125181	HANDYVAC CH14, TROAKAARI (20 KPL)	324619
84010181	KATETERANSATS, STERIL	324627
84010181	UNO CONN CH04-CH18(100/800) INT	324627
40227182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS	324706
40107181	HANDYVAC CH 18 (20/ltk)	324737
84077183	UNOMETER 500 BOLSAS valv. vac 2L. 75uds	324857
40101181	UNO HANDYVAC CH10 BAG W/TAP (1/20) IT	325020
40101181	HANDYVAC CH 10 (20/ltk)	325020
84077182	UNOMETER BAG 2L BO (75/150) IT	325035
84077182	UNO UNOMETER BAG 2L BO (75/150) IT	325035
40103182	UNO HANDYVAC CH12 BAG W/TAP(1/20) INT	325052
40101182	UNO HANDYVAC CH10 BAG W/TAP(1/20) INT	325408
40125183	UNO HANDYVAC CH14 MTL TROC BO(1/20) IT	325409
40223182	HANDYVAC U/BUNDAFTAP M/DRAEN&TROCAR CH14	325410
40105181	HANDYVAC M/DRAEN, CH 14/700 ML POSE	325429

ATTACHMENT 2: EU Product List		
Product Code	Description	Lot Number
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	325429
40121182	UNO HANDYVAC CH 08 MTL TROC BO(1/20) INT	325818
40125181	HandyVac m/draen og trocar CH 14	326205
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	326205
40125181	HANDYVAC CH14, TROAKAARI (20 KPL)	326205
40105181	UNO HANDYVAC CH14 BAG W/TAP (1/20) INT	326619
40123182	UNO HANDYVAC CH10 METAL TROC BO(1/20) IT	326778
40160181	HANDYVAC M/BUNDAFTAPNING CH12	326780
40107181	HANDYVAC M/DRAEN, CH 18/700 ML POSE	327459
40107181	UNO HANDYVAC CH18 BAG W/TAP(1/20) INT	327459
40125181	UNO HANDYVAC CH14 MTL TROC BO(1/20) NEU	327942
40227182	UNO HANDYVAC CH12 MTL TROC WOBO 1X20 AUS	328440