

ASID BONZ GmbH, Hertzstraße 3, 71083 Herrenberg, Germany

Customer name  
 Contact Person  
 Street and house number  
 Postcode and city

ASID BONZ GmbH  
 Hertzstrasse 3  
 71083 Herrenberg, Germany

Phone: +49 7032 9562-0  
 Fax: +49 7032 9562-58  
 Mail: info@asid-bonz.de  
 Web: www.asid-bonz.de

07.04.2020

**Additional warnings for safety information for non-woven disposable gowns (REF 760510 and REF 760520)**

Dear Sirs and Madames,

we ask you to heed and observe the following **Warning**:

Due to identified defects in the sleeve area of the non-woven disposable gowns (REF 760510 and REF 760520), we expressly point out that this product batch cannot guarantee an adequate barrier effect against infectious agents and therefore must be used with special consideration in individual cases.

Please only use in an emergency if no other protective clothing is available. The product is already being revised in order to be able to guarantee the protective effect again in future deliveries.

Productname	REF	PZN	Productcharge	VE	Size (cm, L x W)	Colour
non-woven disposable gowns	760510	11667693	2019-05-31 2019-06-29	10 Stk.	130 x 140	Yellow
non-woven disposable gowns	760520	05886190	2019-07-08 2019-11-04 2017-09-04	10 Stk.	130 x 160	yellow

Please forward this letter to all customers or employees of your facility who use or could use the product. Thank you for your cooperation and understanding.



Marko Rudnicki  
 Geschäftsführer



i.A. Sara Fouladi-Movahed  
 Qualitätsmanagement / Regulatory Affairs